

ADVANTACLEAR®

Surface Disinfectant



Hu-Friedy's AdvantaClear Surface Disinfectant is designed to clean and disinfect dental surfaces while working extremely fast and effective to kill a broad spectrum of disease-causing microorganisms in one minute.*



FAST-ACTING ONE MINUTE CONTACT TIME

A broad-spectrum disinfectant with one minute kill time* – faster kill times mean faster turnover.



SAFE ON SURFACES

Designed to be compatible with a wide range of surfaces including acrylic, aluminum, brass, copper, stainless steel, chrome, plastics and more.



CLEANS & DISINFECTS

Quick and easy way to effectively clean and disinfect surfaces.



FRAGRANCE FREE

Helps eliminate odor-causing bacteria with no harsh fragrance.



LOW ALCOHOL FORMULA

Will not affix blood and soil to surfaces.



Get clear.

How the best perform



* This product kills the following organism in 2 minutes on pre-cleaned, hard, non-porous surfaces at room temperature: Mycobacterium bovis BCG (TB)
See product label for details.

ADVANTACLEAR®

Surface Disinfectant

KILL TIMES vs. LEADING COMPETITOR (in minutes)

		Hu-Friedy AdvantaClear Wipes	Leading Competitor Wipes
Gram-Positive Bacteria	Methicillin-Resistant <i>Staphylococcus aureus</i> (MRSA)	1	2
	<i>Staphylococcus aureus</i>	1	3
	<i>Staphylococcus aureus</i> Rosenbach	1	
	Vancomycin-Resistant <i>Enterococcus faecalis</i> (VRE)	1	2
	Vancomycin Intermediate <i>Staphylococcus aureus</i> (VISA)	1	2
	Methicillin Resistant <i>Staphylococcus epidermidis</i> (MRSE)	1	
Gram-Negative Bacteria	<i>Enterococcus faecium</i>	1	
	<i>Acinetobacter baumannii</i>	1	
	Multi-Drug Resistant (MDR) <i>Acinetobacter baumannii</i>	1	
	<i>Pseudomonas aeruginosa</i> ATCC 15442	1	3
	<i>Pseudomonas aeruginosa</i> ATCC BAA-2108	1	
	<i>Salmonella enterica</i>	1	3
	Carbapenem Resistance <i>Klebsiella pneumoniae</i>	1	
	<i>Klebsiella pneumoniae</i>	1	
	<i>Bordetella pertussis</i>	1	
	Extended-Spectrum β -lactamase <i>Escherichia coli</i> (ESBL)	1	
	<i>Enterobacter aerogenes</i>	1	
	<i>Enterobacter cloacae</i>	1	
	Carbapenem Resistant <i>Escherichia coli</i>	1	
	Enveloped Viruses*	Hepatitis B Virus (HBV)	<1*
Hepatitis C Virus (HCV)		<1*	2
Herpes Simplex Virus Type 2		1	2
Herpes Simplex Virus Type 1		1	2
Human Immunodeficiency Virus (HIV-1) associated with AIDS		<1*	2
Influenza A Virus (H3N2)		1	
Influenza A Virus (H1N1)		1	
Human Coronavirus		1	
Influenza A2 Virus			2
Non-enveloped Viruses*	Rotavirus	1	
	Norovirus	1	
Fungi	<i>Trichophyton interdigitale</i> (ringworm)	1	3
Mycobacteria	<i>Mycobacterium bovis</i> BCG tuberculosis (TB)	2	3

ADVANTACLEAR PRODUCTS

WIPES

Powerful scrubbing wipes that quickly and effectively clean and disinfect while still being safe on surfaces.

- IMS-2160 Canister, 160 Wipes (6 x 6.75 inches)
- IMS-2010 Flat Pack, 10 Wipes (6 x 6.75 inches)
- IMS-2050 Box of Individually Wrapped Wipes, 50 (9 x 8 inches)



READY-TO-USE LIQUID

Convenient, ready-to-use packaging that is easy to use. No mixing or diluting required.

- IMS-2024 Sprayer Bottle, 24 fl. oz. (710 mL)
- IMS-2128 Bottle, 1 Gallon (3.8 L)



Visit us online at Hu-Friedy.com/AdvantaClear

Call 1-800-Hu-Friedy or contact your authorized Hu-Friedy representative for more information.

* HBV, HCV and HIV-1 Precautions: This product kills HBV, HCV and HIV-1 on pre-cleaned environmental surfaces/objects previously soiled with blood/body fluids in health care or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids, and in which the surfaces/objects can be associated with the potential for transmission of HBV, HCV and/or the human immunodeficiency virus Type 1 (HIV-1) associated with AIDS. **Special Instructions for cleaning against HBV, HCV and HIV-1 of surfaces/objects soiled with blood/body fluids: Personal Protection:** Specific barrier protection items to be used when handling items soiled with blood or body fluids are disposable latex gloves, gowns, masks and eye coverings. **Cleaning Procedure:** Blood and other body fluids must be thoroughly cleaned from surfaces/objects before disinfection. **Contact Time:** Allow hard, non-porous surfaces to remain wet for 30 seconds to kill HBV, HCV, and HIV-1 prior to wiping or allow surface to air dry. **Infectious Materials Disposal:** Cleaning materials used that may contain blood/body fluids must be autoclaved and/or disposed of in accordance with local regulations for infectious materials disposal.

