RetraXil[®] blue

KULZER
MITSUI CHEMICALS GROUP

GB	Instructions for use.....	2
FR	Mode d'emploi	4
ES	Instrucciones de uso	6

KULZER
MITSUI CHEMICALS GROUP

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Medical Device – for use by dental health care professionals only.

Intended purpose

RetraXil blue is an astringent retraction paste for the temporary displacement of the marginal gingiva and drying of the gingival sulcus. The paste, which contains aluminium chloride, is inserted into the gingival sulcus. In addition to the displacement effect of the paste, the astringent aluminium chloride ensures a hemostatic effect. Depending on the clinical situation and the practitioner's individual preference, the paste can be used as an alternative or in combination with retraction cords.

Indications

For the temporary displacement of the marginal gingiva and drying of the gingival sulcus, e.g. for

- conventional or digital impressions,
- cementation of temporary and permanent restorations
- and the creation of Class II and V fillings.

Contraindications

The use of this product is contraindicated in case of known or suspected allergies against components of this product. Do not use on patients with periodontal disease, furcation exposure or exposed bone.

Patient target group

This dental accessory is suitable for patients requiring dental treatment for the mentioned indications with consideration of the contraindications. For pregnant and nursing women, the treatment risks must be weighed carefully against the benefits taking into consideration the unborn child or infant.

Intended user

For use by dental health care professionals only.

Clinical benefit

Kulzer medical devices ensure the rehabilitation of oral functions such as chewing, speaking and aesthetics. They stabilize the remaining dentition and/or the alveolar ridge restoratively or prosthetically, or they are accessories in the rehabilitation process.

Potential side effects

Hypersensitivities to the product or its components cannot be excluded in individual cases. Local tissue damage if the application time is too long.

Product features


Primary packaging components and combination products

- Twisting aid (can be used multiple times after sterilization)
- Application cannulas (for single use only – please discard after use)

Composition

Phyllosilicates, water, aluminium chloride hexahydrate, mineral-based oleogel, vinyl siloxane, Pigment Blue 15:3

Product application

Preparation

- Thoroughly rinse the sulcus with water and gently air dry before using the paste. Remove and save the cap (Fig. 1).
- Depending on the clinical situation, a retraction cord can be inserted into the sulcus before the application of the retraction paste.
- Screw on one of the application cannulas supplied (Fig. 2), tighten it using the twisting aid provided and check that it is securely attached (Fig. 3).
- Protect the syringe with application cannula against contamination with a protective sleeve (Fig. 4).
- Apply a small amount of the retraction paste to a mixing pad and discard.

Note: Do not bend the cannula tip.

Application

- Insert the application cannula into the sulcus and apply the retraction paste circularly around the tooth so that excess material wells out of the sulcus (Fig. 5).

Note: Do not use on prepared tooth surfaces – RetraXil blue can impair the adhesive bond of dentine-bonding agents.

- Leave the retraction paste in the sulcus for two minutes (Fig. 6).
– No moisture should be allowed in during this time.

Note: A longer application time can lead to tissue damage.

Rinse / Remove

- Once the application time has ended, remove the retraction paste from the sulcus using an air-water mixture and a suction device; take care not to leave any residue (Fig. 7).
- Check the retraction result and remove any residue from the patient's sulcus/mouth.

Note: If RetraXil blue is not rinsed out of the mouth sufficiently, it can negatively affect the quality of the impression.

Closure

- After application, remove the application cannula along with the protective sleeve using the twisting aid provided (Fig. 8)
- Place the syringe on a hygienically clean surface (Fig. 9)
- Discard the used application cannula along with the protective sleeve (Fig. 10) and reclose the syringe with the supplied cap using clean gloves (Fig. 11)

Note: When using syringes in mouth, use common hygienic protective sleeves to prevent contamination. Do not reuse after contamination.

- The twisting aid supplied must be sterilized after use (Fig. 12). For more information, see **“Cleaning, disinfection and sterilization”**.

Note: Close RetraXil blue carefully after each use.

General information

Please pay attention to the following safety instructions and all processing instructions in other sections.

This product must be used according to the instruction for use only. Any use not complying with these instructions is at the discretion and sole responsibility of the dental professional itself.

Visually inspect for damage before use. Damaged products must not be used. Do not reuse the syringes if there is evidence that the paste has been contaminated with body fluids.

Warnings and safety instructions

- The product may irritate eyes. Avoid contact with eyes. Immediately and thoroughly rinse with plenty of water in case of contact with eyes and contact an ophthalmologist.
- The product may irritate skin. Avoid contact with skin. In case of contact with skin, immediately wash with plenty of water and soap.
- Wear appropriate protective clothing, eye protection and gloves. The patient should use eye protection.
- Do not swallow or ingest this product; this may result in an intolerance.
- Immediately clean paste residues from metal surfaces that have come into contact with the retraction paste, as they can corrode.

Storage conditions

Use within ten weeks of opening the product. Do not store above 25°C (77°F). Store in a dry place. Do not use the material after the expiry date.

Disposal information

Recommendation: Dispose of in accordance with official regulations. Do not dispose of contents or partially emptied packaging in the household waste or allow them to enter the sewage system.

Cleaning and disinfection

WARNING: Please observe the country-specific regulations on hygiene standards and on the cleaning/disinfection of medical devices in dental practices. Wear suitable gloves. The twisting aids are supplied unsterilized. The twisting aid must be sterilized prior to first use and cleaned, disinfected and sterilized after each use.

Processing restrictions

Oxidizing disinfectant solutions are not suitable. Strong acids and strong bases can damage the plastic.

TWISTING AID instructions

In treatment room

Unscrew the twisting aid from the syringe and discard the cannula along with the protective cover.

Retention and transport

For transport outside of treatment room use a closed container. It is recommended that the twisting aid be reprocessed as soon as possible, but no more than two hours after use.

Preparing for decontamination

Wear personal protective equipment (sturdy gloves, water-repellent protective gown, face mask or safety goggles and mask).

A) Cleaning and disinfection: by hand

Equipment: non-protein-binding locally approved virucidal disinfectant with cleaning action based on quaternary compounds (e.g. CaviCide™), lint-free cloths, plastic brush, if possible fully demineralized water (free from facultative pathogenic microorganisms), tub for cleaning agents and disinfectants.

- 1) Prepare the cleaning solution according to the manufacturer's instructions.
- 2) Immerse the twisting aid completely in the cleaning solution.
- 3) While immersed, brush the hard-to-reach areas of the immersed twisting aid with a soft brush. Pay attention to the critical hard-to-reach areas and cavities where it is not possible to visually assess the cleaning effect.
- 4) Immerse in the cleaning solution according to the manufacturer's instructions.
- 5) Remove the twisting aid from the cleaning solution and rinse thoroughly (e.g. 30 seconds) under running water.
- 6) Check for cleanliness; if contamination is still visible, repeat the above steps.
- 7) Prepare the disinfectant according to the manufacturer's instructions.
- 8) Immerse the twisting aid completely in the disinfectant.
- 9) Immerse in the disinfectant according to the manufacturer's instructions (e.g. 5 minutes).
- 10) Remove the twisting aid from the disinfectant and rinse thoroughly (e.g. 30 seconds) with demineralized water.

11) Dry with a lint-free disposable cloth and/or with medical compressed air.
Only validated processes may be used for the cleaning/disinfection of twisting aids.

B)Cleaning and disinfection: Automated WD [washer disinfecto

Equipment: cleaning and disinfection device with thermal program (temperature 90°C/194°F to 95°C/203°F),
Cleaning agent: mildly alkaline cleaning agent.

- 1) Brush the hard-to-reach areas of the twisting aid for one minute under running potable water with a soft brush. Pay attention to the critical hard-to-reach areas and cavities where it is not possible to assess the cleaning effect visually.
- 2) Place all twisting aids in a suitable mesh tray or place them on the load carrier so that all internal and external surfaces of the instruments are cleaned and disinfected.

Program step	Water	Amounts required	Time	Temperature
Pre-wash	CW (cold water)		5 min	
Cleaning solution dosing		Acc. to the manufacturer's specifications		Acc. to the manufacturer's specifications
Cleaning	DW (fully demineralized water)		10 min	55°C (131°F)
Rinsing	DW (fully demineralized water)		2 min	
Disinfection	DW (fully demineralized water)		3 min	A0 value > 3000 ¹ (e.g. 90°C / 194°F, 5 min)
Drying			15 min	Up to 120°C (248°F)

¹ The authorities may issue other parameters for disinfection performance in their area of responsibility.

- 3) Remove all twisting aids at the end of the program.
- 4) Check that the load is dry and, if necessary, dry with medical compressed air or with a lint-free cloth.
- 5) After removal from the washer disinfecto

Maintenance, inspection and testing

All instruments are to be checked visually for cleanliness, integrity and functionality, if necessary, using an illuminated magnifying glass (3-6 dioptre).

All twisting aids must be checked for damage and wear. Damaged products must not be reused and must be disposed of.

C) Sterilization

Device: Sterilizer or small steam sterilizer.

Packaging:

- Equipment: Foil-paper packaging, sealing device
A suitable procedure (sterile barrier system) must be used for packaging the twisting aids.
- Individual: A sterile barrier system (foil-paper packaging PeelVue + sterilization pouch) should be used, one that is intended by the manufacturer for steam sterilization. The packaging must be large enough so that the sealing seam is not under tension. (During the validation, double packaging is carried out in foil-paper packaging).
- Note: After the heat sealing process, the sealing seam must be checked visually for possible defects. In the event of errors, the packaging must be opened and the instrument repackaged and sealed.
- Sets: Sort twisting aids into the appropriate sieves so that they cannot be mechanically damaged during transport. Sterile barrier systems in accordance with ISO 11607 (e.g. sterilization fleece) (the instruments can be single packed or double packed) or sterile goods containers are suitable for packaging the sieves. The weight of the sterile goods container must not exceed 10 kg.

Procedure: Steam sterilization with fractional pre-vacuum, 134°C (273°F), holding time at least 3 min or 132°C (270°F) at least 4 min.
Drying time: 15 minutes.

- 1) Place the packaged twisting aid in the sterilization chamber.
- 2) Start the program.
- 3) Remove the products at the end of the program and allow them to cool.

Check the packaging for any damage or moisture penetration. Compromised packaging is to be assessed as non-sterile. The instruments must be repackaged and sterilized.

Storage: Storage (protect from dust and moisture) and storage duration at the user's premises according to the specifications.

Only validated processes may be used for the sterilization of twisting aids.

Complaint notice

Please quote batch number and article number or UDI in all correspondence about the product.
Article number REF and batch number LOT or UDI: See note on syringe/packaging.

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 +49 800 4372522

Made in Germany

Distributed in USA/Canada
exclusively by:

Kulzer, LLC
4315 South Lafayette Blvd.
South Bend, IN 46614-2517
1-800-431-1785

Caution:
Federal law restricts this device
to sale by or on the order of a dentist.

Explanation of symbols on labelling

	Medical device
	Reorder no.
	Manufacturer
	Batch code
	Use-by date
	Date of manufacture
	Consult instructions for use
	Storage temperature
	Keep dry
	Do not re-use
	Unique Device Identification

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Kulzer GmbH