

LASCOD

dental clinic use



Ghenesyl

addition curing polyvinylsiloxane silicone

Polyvinylsiloxane addition curing silicone is a versatile and efficient product guaranteeing maximum precision in humid conditions. Ghenesyl, now with its renewed formula, has different grades of hardness (Hard/Soft) and consistency (Superlight/Light/Regular/Heavy) which attain the best results adapting to every impression technique.

A long working time increases the product workability. At the same time a short setting time in mouth is guarantee for the best patient's comfort.

new

FastSet

Ghenesyl is now available in two different setting releases: normal and the all-new fast version which is faster in the working time and setting time in mouth.

mint flavour

25µm

Ghenesyl Putty, Normal Set and Fast Set are available now with improved workability and a slight mint flavour.

PUTTY



New formula!

5µm

thixotropicity,
accuracy
& hydrophilicity




BODY



Ghenesyl, now with its renewed formula, is highly thixotropic and hydrophilic. The material stays in place after dispensing and its flows well under pressure. Ghenesyl easily adapts to the anatomical shape of dental structures and soft tissues, thus reproducing perfectly integral margins.

<10°



The contact angle measurement (wettability index) is less than 10°.



Putty
Mint Flavour



GHENESYL KIT

Ghenesyl is available in KIT PACKS containing 2x150ml putty base and catalyst and a cartridge of 50ml body.

	Working time including mixing time (23°)		Minimum time in mouth (37°)		Recovery from deformation %	Dimensional stability	PUTTY/ WASH	SANDWICH	MONOPHASE
	NORMAL SET	FAST SET	NORMAL SET	FAST SET					
PUTTY SOFT	2' 30"	1' 15"	2' 15"	1' 30"	99.9	7 days			
PUTTY HARD	2' 30"	1' 15"	2' 15"	1' 30"	99.9	7 days			
SUPERLIGHT	2'	1' 20"	1' 30"	60"	99.9	7 days			
LIGHT	2'	1' 20"	1' 30"	60"	99.9	7 days			
REGULAR	2'	1' 20"	1' 30"	60"	99.9	7 days			
HEAVY	2'	1' 20"	1' 30"	60"	99.9	7 days			
MONO	2'	-	1' 30"	-	99.9	7 days			

Ghenesyl mono

monophase addition curing
polyvinylsiloxane silicone

Ghenesyl Mono, developed by LASCOD, is a new impression material with a medium consistency, ideal for monophase technique. Recommended for crowns and bridges impression taking, implantoprosthesis, inlay, onlay, partial and total edentulism. Thanks to its characteristic viscosity it can be used directly on the surfaces to be reproduced and on impression tray. It's syringable through the included mixing tips or through the specific syringe for elastomers. Given its optimal flexibility Ghenesyl Mono can be easily removed from the mouth without structural alteration.

Ghenesyl Mono can be used also with Super Light Body / Light Body / Regular Body silicones as support material inside individual impression tray or standard stainless steel impression tray.

- Ideal viscosity for monophase technique
- High details definition (5 μm)
- Maximum thixotropicity and hydrophilicity
- Fast setting time in mouth
- High resistance to tear
- High elastic recovery
- High dimensional stability over time
- Slight mint flavour
- Resistant to aseptic cleaners



5 μm

Mint
Flavour



Oklurest

silicone for
bite registration

Oklurest is an addition silicone (polyvinylsiloxane) made for orthodontic bite registration, registration keys for gnathologic registrations, inter-maxillary registration keys for centered positions, eruptions and ectopic eruptions, registration for cephalometric analysis, which requires a CAD scanning (optical, laser and tactile).

Oklurest characteristics grant a working base which allow the technicians to obtain the maximum precision.

The high dimensional stability over time ensures optimal conservation even in the presence of temperature changes.

imperceptible viscosity

Oklurest is easy-to-use and offers a series of advantages which allow the maximum control during positioning.

The viscosity is imperceptible for patient.



The useful hardness degree (95 Shore A) achieved in short time, allows to work on masses widely stable facilitating the finishing stages by cutters or burs.

Working time (23°C)	45"
Minimum time in mouth (37°C)	60"
Hardness Shore A	95
Linear dimensional change (after 24 hours)	-0,02%



- Easily and safely workable thanks to double cartridge system 1:1
- Maximum control during positioning
- Thixotropic
- Imperceptible viscosity for patient
- Exact reproduction of occlusal details (5µm)
- Short setting time
- Minimum dimensional variation
- Undeformable / stable over time
- Easy insertion and removal from the model
- Easy finishing by burs

maximum precision



5µm

95 Shore A

**CAD
SYSTEM**

Silaxil

condensation curing silicone

The last technological development of C-silicone at Lascod providing the highest accuracy in details survey. Adaptable to every impression technique. Its high dimensional stability allows preserving the impression over time without any dimensional alteration of the models.

40µm

PUTTY

5µm

LIGHT BODY

5µm

REGULAR BODY



Silaxil is now available also in the convenient Silaxil Box and the easy-to-use Silaxil Kit (refer to the last page' table for the item codes).





Its consistency ease mixing, making it free from any porosity, it also allows the detail's reproduction without compressing the mouth tissues.

accuracy

fast setting time




The long working time allows to perform the various stages of impression taking, without the material pre-hardens hastily.

Minimum time in mouth!

Elastic returns as well as compressive strength are well beyond the standards of the products currently on the market.

flavours & colours



The selection of delicate flavours, mint for the Putty and orange for the Light Body and the Regular Body, allows not to alter the degree of the patient's saliva and greatly reduces the emetic effect.

The purple colour (Light Body) and the green colour (Regular Body) allow to check the perfectly homogeneous mixing of base paste and the catalyst, creating later a great contrast between the first and the second impression.
















Moreover, thanks to the blocks mixture driven, the phases of the dosage and the subsequent mixing between Light Body/ Regular Body and the Energyl catalyst are extremely simplified, avoiding a useless waste of material.



	PUTTY	LIGHT BODY	REGULAR BODY
Mixing time (23°C)	30''	30''	30''
Working time (23°C) including mixing time	1' 20''	2'	2'
Minimum time in mouth (37°C)	3'	3' 30''	3' 30''
Recovery from deformation	99.75%	99.9%	99.9%
Strain in compression	3.5%	4.5%	4.5%
Hardness Shore A	50 SHA	40 SHA	40 SHA

packaging

*** ALSO AVAILABLE
IN ECONOMY PACK
4 x 50 ml cartridges
without mixing tips**

GHENESYL		NORMAL	FAST	
	PUTTY HARD	GNS020	GNS520	2 x 300 ml base + catalyst
	PUTTY SOFT	GNS010	GNS510	2 x 300 ml base + catalyst
	HEAVY BODY*	GNS070	GNS570	2 x 50 ml cartridges + 12 mixing tips
	SUPERLIGHT BODY*	GNS040	GNS540	2 x 50 ml cartridges + 12 mixing tips
	LIGHT BODY*	GNS050	GNS550	2 x 50 ml cartridges + 12 mixing tips
	REGULAR BODY*	GNS060	GNS560	2 x 50 ml cartridges + 12 mixing tips
	MONO* (50 ml)	GNS400	-	2 x 50 ml cartridges + 6 mixing tips
	MONO (380 ml)	GNS480	-	2 x 380 ml cartridges + 2 mixing tips
	KIT PUTTY HARD + LIGHT BODY	GNS080	GNS580	
	KIT PUTTY HARD + SUPERLIGHT BODY	GNS081	GNS581	2 x 150 ml (putty) + 1 x 50 ml cartridge (body) + 6 mixing tips + 6 oral tips
	KIT PUTTY SOFT + LIGHT BODY	GNS082	GNS582	
	KIT PUTTY SOFT + SUPERLIGHT BODY	GNS083	GNS583	
OKLUREST				
	<i>Silicone for bite registration *</i>	OKR095	-	2 x 50 ml cartridges + 12 mixing tips
SILAXIL				
	PUTTY	SLP010		1 x 900 ml
	LIGHT BODY	SLP050		1 x 140 ml
	REGULAR BODY	SLP060		1 x 140 ml
	ENERSYL Catalyst	ENS060		1 x 60 ml
	SILAXIL BOX 1 x 900 ml (putty) + 1 x 140 ml (light body) + 1 x 60 ml (catalyst Energyl) + 10 pcs of mixing guide	SLP301		
	SILAXIL KIT 1 x 200 ml (putty) + 1 x 40 ml (light body) + 1 x 15 ml (catalyst Energyl) + 3 pcs of mixing guide	SLP201		
		SLP150 Mixing guide (booklet of 30 pcs)		
SILICONE ACCESSORIES				
				
		Dispenser (50ml) GNS100		48 mixing tips (50ml) GNS200 - superlight/light GNS250 - regular/mono GNS270 - heavy - OkluREST
		60 oral tips GNS300 - superlight/light GNS350 - regular/heavy		24 mixing tips (380ml) GNS280 - mono

The photographic images here reproduced are purely indicative and are not necessarily identical to the actual products.

LASCOD

LASCOD S.p.A. - Via L. Longo, 18 - 50019 Sesto Fiorentino, Florence, Italy - tel.: +39 055 4215768 - e-mail: lascod.italy@lascod.it
www.lascod.com

For further
information
activate the
QR Code reader
on your smartphone.

