



Limit uncertainty with quality instrumentation
from the brand you trust.

Xcision® Crown Spreader Forcep

A member of our Xcision® line, this stainless steel crown spreader quickly splits a crown resulting in its rapid removal.

Simply score the crown with a bur creating a slot wide enough for the crown spreader beaks. Place the crown spreader into the slot and squeeze the handle allowing the beaks to separate and split the crown.



DEFXCCS

Wynman Crown Gripper

The tungsten carbide inserts on this instrument will help grip the temporary crown securely without slipping. This secure grip aids in the removal of the temporary crown while still intact.



72-10TC

Morrell Crown Remover



Crown Keys



Xcision® Crown Spreading Forcep

Reference	Description
DEFXCCS	Xcision® Crown Spreading Forcep

Morrell Crown Remover

Reference	Description
72-2	Morrell Crown Remover with all 3 tips
72-4A	Point A only
72-4B	Point B only
72-4C	Point C only

Wynman Tungsten Carbide Crown Gripper

Reference	Description
72-10TC	Wynman Tungsten Carbide Crown Gripper

Crown Keys

Reference	Description
KKL	Crown Key Long
KK	Crown Key Straight
KK90	Crown Key 90°
KK45	Crown Key 45°

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instructions for use.

For more information or to place an order, please contact:

United States, Canada, Asia, Pacific, Latin America
 USA 800-654-2873 ■ 888-980-7742 fax
 International +1 609-936-5400 ■ +1 609-750-4259 fax
integralife.com/contact

Integra York PA, Inc.
 589 Davies Drive
 York, PA 17402 ■ USA