

Alcohol Torch

Economical alcohol torch that produces a needle-point flame



The classic Buffalo Dental Alcohol Torch produces a needle-point flame with only a gentle squeeze of the plastic bottle, allowing for precision control.

The Alcohol Torch now features an optional 3-point base, which reduces the chance of knocking the bottle over.

Features and Benefits:

- 5 1/4" high torch uses denatured alcohol fuel for cost savings and convenience
- Base of bottle is weighted with marbles to further improve safety
- Ideal for use in applying rolled edges to vacuum formed adaptations (flame edge and rub with thumb)
- Versatile torch can also be used as an alcohol lamp
- 3-point base for added stability
- The base has an easy slide fit and can remain with the torch as you use it or on the bench as a receptacle
- Includes plastic cap to help keep fuel from evaporating when torch is not in use (for use after torch is extinguished)



Ordering Information:

| Item# | Description |
|--------------|--|
| 82440 | Plastic Alcohol Torch w/ 7" Wick |
| 82410 | Plastic Alcohol Torch w/ 7" Wick and Base |
| 57560-1 | Replacement Wick, 7" long, Pkg. of 1 |
| 57560-1D | Replacement Wicks, 7" long, Pkg. of 12 |
| 82400 | Fluorescent Base Only |
| | Torchcover Replacement Clear Cap for Torch |

Buffalo Dental Manufacturing Co., Inc.

159 Lafayette Dr., PO Box 678, Syosset, NY 11791-0678 | Tel: 516.496.7200 or 800.828.0203

Fax: 516.496.7751 | Email: sales@mailbdm.com | www.buffalodental.com