


SciCan
HYDRIM

Wash Test Indicators

Monitoring Instrument Washers



Think of it as a routine check up for your washer

The **HYDRIM Wash Test Indicators (WTIs)** have been designed to provide two-dimensional validation of HYDRIM washing efficacy. It is a simple, easy-to-use indicator that confirms instruments are effectively and consistently cleaned with HIP Cleaning Solution.



Safe

Synthetic test soil (non-toxic indicator) designed to monitor the efficacy of HYDRIM washers.



Reliable

Designed to provide verification of a consistently effective cleaning process in HYDRIM washers.



Smart Design

Indicators are folded “back-to-back” and placed inside the holder, monitoring spray performance from two directions.




Easy

WTI wash off characteristics clearly indicate whether the wash cycle has passed or failed; clean indicator at the end of the cycle indicates an effective cleaning process.

 **COLTENE**

The only validated and effective cleaner in the market for use in HYDRIM units.

HIP Cleaning Solution is designed to remove any type of protein/debris from instruments. HIP is automatically added at various stages in the cycle to optimize the cleaning and protection of instruments. To maximize HYDRIM cleaning process, it is essential to use the certified and validated HIP Cleaning Solution that works for your HYDRIM unit during the wash cycle.



HIP Cleaning Solution for HYDRIM L110W

Designed for HYDRIM L110W, HIP Cleaning Solution box fits inside L110W, instead of taking up cabinet space. Strong, built-in tube connector prevent leaks.

HIP Ultra Cleaning Solution for HYDRIM C61W

Designed for HYDRIM C61W, HIP Ultra Cleaning solution pouches fit inside the C61W, instead of taking up cabinet space. Strong, built-in tube connector prevent leaks.

Part Numbers

PART NUMBER	DESCRIPTION
WTI-SK	Wash Test Starter Kit (One Holder and 32 Indicators)
WTI-RK	Wash Test Refill Kit (96 Indicators)
CS-HIPL	HIP Cleaning Solution for L110W G4 (2x3.8L)
CS-HIPC-U	HIP Ultra Solution for C61W/WD G4 (8x750mL)
C61W-D01	HYDRIM C61W G4
L110W-D01-G4	HYDRIM L110W G4

Manufactured for:
SciCan Ltd.
 1440 Don Mills Rd.
 Toronto, ON
 M3B 3P9 Canada

Distributed by:
SciCan Ltd.
 1440 Don Mills Rd.
 Toronto, ON
 M3B 3P9 Canada

Distributed by:
Coltene/Whaledent Inc.
 235 Ascot Pkwy
 Cuyahoga Falls, OH
 44223 USA

www.scican.com

