Dux Dental Material Safety Data Sheet

 Dux Dental
 Phone: 800.833.8267

 600 E. Hueneme Road
 Fax: 805.488.2266

 Oxnard, CA, US
 Fax: 805.488.2266

 Product
 Gingibraid, Gingiknit, Gingiyarn, Unibraid, Uniknit & ShortcutAlum

 Other Name
 Aluminum Potassium Sulfate Gingival Retraction Materials

 REF
 13076, 13077, 13078, 13210, 13211, 13212, 13213, 13360, 13361, 13362, 13363, 13365, 13366, 13367, 13368, 13508, 13509, 13510, 13511, 13512, 13513, 13608, 13609, 13610, 13611, 13612, 13613, 13730, 13731 & 13732

MSDS Number MS113

Hazardous Ingredients				
Principal Hazardous Components None		pel / tlv	CAS Number NA	Percentage NA
		NA		
Physical and Chemical Charat				
Appearance	Small cords of varying diameters and con	nstructions impregnated v	v/ Alum.	
Boiling Point	NA			
Melting Point	NA			
Vapor Pressure	NA			
Evaporation Rate	NA			
Odor	None			
Vapor Density	NA			
Specific Gravity	NA			
Water Solubility	Yes (Alum)			
Fire & Explosion Data				
Explosion	Not considered to be an explosion hazard	1.		
Fire	Not considered to be a fire hazard.			
Fire Ext. Type	А			
Flashpoint	NA			
Special Fire Procedures	None			
Unusual Hazards	None Known			
Reactivity Data				
Stability	Stable			
Incompatibilities	Sl. corrosive to metals.			
Hzd Decomp Products	Oxides of Sulfur			
Hazardous Polymerization	Will not occur.			

Dux Dental Material Safety Data Sheet

Dux Dental 600 E. Hueneme Road Oxnard, CA, USA

 Product
 Gingibraid, Gingiknit, Gingiyarn, Unibraid, Uniknit & ShortcutAlum

 Other Name
 Aluminum Potassium Sulfate Gingival Retraction Materials

 REF
 13076, 13077, 13078, 13210, 13211, 13212, 13213, 13360, 13361, 13362, 13363, 13365, 13366, 13367, 13368, 13508, 13509, 13510, 13511, 13512, 13513, 13608, 13609, 13610, 13611, 13612, 13613, 13730, 13731 & 13732

MSDS Number MS113

Health Hazard Information				
Ingestion	NA			
Inhalation	Dust may cause irritation.			
Skin/Eye Contact	Dust may cause irritation. Flush area w/ water if exposed.			
Toxicity Data	May cause reversible or irreversible changes to exposed tissue, not permanent injury or death; can cause considerable discomfort. Use in accordance with manufacturers instructions and sound professional			
Acute Hazards	None			
Chronic Hazards	None			
Carcinogenicity	NA, no product ingredient listed by NTP, IARC, or OSHA.			
Precuations for Safe Handling and Use				
Cleanup	NA			
Disposal	Dispose of in accordance with local, state, and federal regulations.			
Precautions	As part of good industrial, personal hygiene, and safety practice avoid all unecessary exposure to the substance and ensure prompt removal from skin, eyes, and clothing.			
Occupational Control Measure	28			
Eye Protection	NA			
Gloves	NA			
Dust/Particulate Mask	NA			
Other	NA			
Additional Information				
DOT Hazard Class	Not Regulated			
DOT Small Qty Exception	NA			
DOT Other Exception	NA			
DOT Exception Permitted	NA			
Allergenic Ingredients	This product does not contain Latex, Peanut Products, Gluten, or Red Dye #3.			

The information contained herein is provided in good faith and is believed to be correct as of the date hereof. However, the company makes no representation as to the comprehensiveness or accuracy of the information. It is expected that individuals receiving the information will exercise their independant judgement in determining its appropriateness for a particular purpose. Accordingly, we will not be responsible for damages of any kind resulting from the use of or reliance upon such information. No representations or warranties, either express or implied, of merchantability, fitness for a particular purpose or of any other nature are made hereunder with respect to the information set forth herein or to the product to which the information refers.