

# Innovation You Can See. Performance You Can Trust.

## IntelliTek® Technology

Unlike traditional light-cured composites that shrink toward the curing light, BulkEZ Plus utilizes dual-cure technology to control and direct shrinkage, helping reduce gap formation and microleakage. Its patented initiator system delivers reliable bond strength without the need for dual-cure activators.



### Unlimited Depth of Cure

Restore deep preparations confidently without layering restrictions, saving valuable chair time.



### 80 nm Spherical Zirconia Fillers

Delivering enhanced polishability, wear resistance, and improved esthetics compared with the original formulation.



### One-Step Simplicity

Injectable handling, rapid dual-cure, exceptional esthetics, and high strength deliver true quadrant dentistry without compromise. Compatible with universal bonding agents, BulkEZ Plus is a true one-step bulk-fill composite—no dual-cure activator or capping layer required.



### Built for More Than Posterior Restorations


BulkEZ Plus is indicated for:

- Class I-VI restorations
- Crown fabrication
- Core build-ups
- Post cementation



### One Composite. Multiple Indications.

From routine posterior composites to complex restorative procedures, BulkEZ Plus delivers the versatility clinicians need with the simplicity they want.



Advanced Nanofillers and uniform clusters make high loading and adaptability possible




Courtesy of Dr. Troy Schmeding.

# Enhanced esthetics. Unlimited cure depth. Proven IntelliTek® Technology.

Restore with confidence, efficiency, and predictable results—case after case. BulkEZ Plus delivers the performance clinicians trust.

Available through distributors nationwide and from Zest Dental Solutions.

## bulkeZ PLUS® Here are the Facts:



\* Source: Evaluation of Various Properties of "Bulk EZ Prototype" Oregon Health & Science University

**ZD ZEST DENTAL  
SOLUTIONS®**

Zest Dental Solutions  
2875 Loker Ave E | Carlsbad, CA 92010  
1-800-262-2310 • Zestdent.com

L8166 REV A 07/2026



Scan To Learn More

**ZD ZEST DENTAL  
SOLUTIONS®**



**bulkeZ PLUS®**

## In a Class of Its Own... the Bulk Fill Composite that Delivers!



## Meet the BulkEZ Plus® Family

This complete bulk-fill composite solution delivers:

Unlimited depth of cure

Enhanced esthetics

Expanded shade selection

Designed to simplify restorative workflows while delivering the strength, polishability, and reliability clinicians expect from BulkEZ Plus®.

## One System. More Possibilities.

BulkEZ Plus combines the simplicity of a dual-cure bulk-fill composite with advanced esthetics and long-term durability. Powered by IntelliTek® Technology and enhanced with 80 nm spherical zirconia fillers, clinicians can restore with confidence, whether treating a routine posterior or a complex restorative case.



## The Clinician's First Choice for Performance:

- Unlimited depth of cure
- Directed shrinkage technology achieving ZERO microleakage\*
- Exceptional polishability and wear resistance
- Injectable composite offering the adaptability of a flowable with the filler loading of a packable composite
- Broad indications from direct restorations to core build-ups and post cementation

\*Source: Microleakage under Class II restorations restored with various bulk-fill materials. University of Alabama at Birmingham

## One Family. Every Shade. Every Procedure.

BulkEZ Plus now offers a complete shade lineup to support today's esthetic demands while maintaining the simplicity of a single restorative system.

### Available Shades:



## Exceptional Esthetics. Proven Strength.

Advanced spherical zirconia filler technology delivers life-like esthetics, exceptional polish retention, and seamless shade integration—while maintaining the strength and durability clinicians expect for long-lasting restorations.



Flexural Strength: 146 MPa



Compressive Strength: 286 MPa

## Precision Placement. Your Way.

BulkEZ Plus is available with multiple cannula options to help clinicians achieve optimal access, control, and material placement based on the procedure and personal preference.


### More Material. More Value.

Each 6g syringe delivers 50% more product than the industry-standard 4g syringe.



| Cannula  | Features & Benefits                                                  | Best For                            |
|----------|----------------------------------------------------------------------|-------------------------------------|
| 17 gauge | Widest - Fast placement for large restorations.                      | Large restorations & core build-ups |
| 19 gauge | Balanced - The perfect combination of flow and control.              | Everyday restorative procedures     |
| 20 gauge | Precise - Controlled placement in deep or conservative preparations. | Deep access & precise placement     |

## Application Guide



Always review the IFU before using the products and follow the steps as instructed in the IFU.