



# ENSURE™

Monitoring for Steam Sterilization

## Think of it as a routine check up for your sterilizer.

The **Ensure**™ line of Infection Control monitoring products ensures your equipment is performing at its peak. It includes a range of products that are simple, easy to use and help to ensure your patients are being kept safe with every cycle.



### Safe

Designed to help keep your patients safe with every cycle.



### Easy

Provide the highest level of sterility assurance with easy-to-read indicators, which clearly identify if an instrument load met the expected parameters for sterilization.



### Capacity

The 10-well incubator allows you to monitor up to 5 sterilizers at one time, providing you with easy-to-read results that show whether sterility has been achieved or not.



### Reliable

Meets the most stringent requirements in today's changing regulatory environment.

# Ensure™ Sterilization Monitoring System

## Ensure™ Accessories



### Ensure™ Chemical Indicators and Integrators

Ensure™ Type 4 and Type 5 indicators are fast and economical solutions for load monitoring. Type 5 Integrators are the most accurate of all chemical indicators and contain a chemical ink which reacts to all three sterilization parameters (time, temperature and pressure).

### Ensure™ Sterilization Pouches

Ensure™ Sterilization Pouches have built-in multi-parameter internal/external process indicators that change color when all three criteria for sterilization have been met - time, temperature and pressure. Once sterilized, contents will remain sterile for up to 3 years provided the integrity of the pouch has not been compromised.

### Ensure™ Biological Indicators (BI)

Ensure™ BIs are reliable and self-contained indicators for results in as little as 24 hours. BIs are the only method to ensure effective performance of your sterilizer. The Ensure™ BIs are available in a starter kit system that consist of 25 BIs, the unique incubator and an easy to fill-in log book.

### Ensure™ Incubator

The Ensure™ Incubator is a unique 10-well incubator. With an integrated vial crusher and built-in timer for up to five different programs, the SciCan Ensure™ BI system leaves no room for error.

### Part Numbers

Part Number	Description
SCI134	Type 4 Chemical Indicators (250/box)
SCI-T5	Type 5 Chemical Integrators (100/box)
SCI-BI100	Biological Indicators (100/box)
SCI-BISK-120	BI Starter Kit (25 BI's, incubator, log book)
SCIBINC-120	Ensure™ Incubator
SCI-SCP	Sterilization Pouches (200/box)

Distributed by:  
**SciCan Inc.**  
701 Technology Dr. Canonsburg,  
PA 15317 USA  
Tel: 724-820-1600  
Fax: 724-820-1479  
Toll free: 1-800-572-1211  
uscustomerservice@scican.com

Manufactured for:  
**SciCan Ltd.**  
1440 Don Mills Road  
Toronto ON M3B 3P9  
CANADA  
Tel.: 416-445-1600  
Toll free: 1-800-667-7733  
Fax: 416-445-2727