



# STATMATIC

Handpiece Maintenance Unit

## Ultimate Care System for handpieces and their chuck mechanisms.

**STATMATIC** is the automatic care and maintenance system that properly cleans, lubricates, and purges up to 3 handpieces prior to sterilization.

The STATMATIC also incorporates a unique chuck cleaner to dislodge debris inside the chuck.

You and your staff will profit from greater efficiency and fewer repairs of your handpieces.



### Extends the Life of your Handpieces

It properly cleans, lubricates, and purges your handpieces to increase their working lifetime.



### Optimizes Staff Productivity

Its automated process optimizes working procedures in your practice, allowing staff to focus on alternative tasks.



### Time-Saving

Up to 3 handpieces can be cleaned, lubricated, and purged in less than a minute.



Compressed air connection required

## Technical Specifications

|                               |   |
|-------------------------------|---|
| <b>Dimensions (L x W x H)</b> | 28.5 x 19 x 40 cm / 11.4" x 7.6" x 16"  |
| <b>Weight</b>                 | 8.6 kg / 18.9 lbs   |
| <b>Electrical Ratings</b>     | Rated voltage: 100-230 V<br>Rated frequency: 50/60 Hz<br>Rated power: max. 30 VA  |
| <b>Protection Class</b>       | 1   |
| <b>Ambient Temperature</b>    | 5-40° C / 41-104 °F   |
| <b>Sound Level</b>            | 60 (65), dB (A)   |
| <b>Compressed Air</b>         | 40 NL/min   |
| <b>Air Consumption</b>        | max. 85 %   |
| <b>Compressed-Air Quality</b> | Residual oil content: max. 0,1 mg/m <sup>3</sup><br>Residual dust: max 1µ; max. 1 mg/m <sup>3</sup><br>Residual water: max. 0.1 g/m <sup>3</sup> at -40°C |

# Automatic Care & Maintenance

for up to 3 instruments  
plus their chuck mechanisms

**Automated, reproducible quality of care in less than a minute.**



Dedicated chuck care port for maintaining the handpiece chuck



Easy spray can insertion with unique safety valve



Absorption pad for excess lubricant and debris

Manufactured for:  
**SciCan Ltd.**  
1440 Don Mills Road,  
Toronto, ON M3B 3P9  
Canada  
T 416-445-1600  
F 416-445-2727  
TF 1-800-667-7733

STATMATIC is a trademark of SciCan Ltd.

Distributed by:  
**COLTENE Inc.**  
235 Ascot Parkway  
Cuyahoga Falls, OH 44223  
USA  
T 330-916-8800  
info@coltene.com

For more information



Call Us  
**+1 416-446-1600**  
Toll Free: **+1 800-667-7733**



Email Us  
**customerservice@scican.com**

[www.scican.com](http://www.scican.com)

**SciCan** | Your Infection Control Specialist