

SEE YOUR WORLD IN A DIFFERENT LIGHT



DARK MODE:  ON

Stars shine brightest when it's dark Be an Absolute Star.

Your confidence in your glove's performance to comfortably keep you and your patient safe should be Absolute. Introducing the Aurelia Absolute with a special soft formulation 3.2mil construction and our Aurelia Premium Protection suite of quality and performance enhancing features to give you a new sense of Absolute precision and control as you go through your day.



AURELIA® GLOVES
For samples visit aureliagloves.com



Absolute

Product Description

Aurelia Absolute Black Nitrile Gloves features a new soft nitrile formulation for all day comfort and ergonomics, finished in a high contrast satin black finish and featuring a full finger texture for enhanced grip and precision control. Aurelia Absolute is ideal for light to medium applications requiring high finger sensitivity, tactility, flexibility and optimal barrier performance. www.aureliagloves.com



Quality Signs & Certificates

89/686/EEC	EN 420 : 2003 A1 : 2009
93/42/EWG	EN 374 2-3
ISO9001:2008	EN 388
ISO 13485	EN 1186 : Food Contact Approved (SGS)
CMDCAS 13485 (SGS)	ASTM D 5250
EN 455 1-4	AQL 1.5

Aurelia® Absolute® Specifications

TENSILE TEST RESULTS (unaged)

	Tensile Strength (MPa)	Elongation at break (%)
Extra Small	min 14	min 500
Small	min 14	min 500
Medium	min 14	min 500
Large	min 14	min 500
Extra Large	min 14	min 500
ASTM D6319-00	14 MPa minimum	500% minimum

AGED TENSILE TEST RESULTS (7 Days at 70° C)

	Tensile Strength (MPa)	Elongation at break (%)
Extra Small	min 14	min 400
Small	min 14	min 400
Medium	min 14	min 400
Large	min 14	min 400
Extra Large	min 14	min 400
ASTM D6319-00	14 MPa minimum	400% minimum

DIMENSIONS

	Length (inches/mm)	Thickness (Single-wall)	
		Palm (mils/mm)	Finger (mils/mm)
Extra Small Aurelia Absolute ASTM D6319-00	9.45 / 240 8.7 / 220 minimum	3.2 / 0.08 1.96 / 0.05 minimum	3.94 / 0.10 1.96 / 0.05 minimum
Small Aurelia Absolute ASTM D6319-00	9.45 / 240 8.7 / 220 minimum	3.2 / 0.08 1.96 / 0.05 minimum	3.94 / 0.10 1.96 / 0.05 minimum
Medium Aurelia Absolute ASTM D6319-00	9.45 / 240 9.1 / 230 minimum	3.2 / 0.08 1.96 / 0.05 minimum	3.94 / 0.10 1.96 / 0.05 minimum
Large Aurelia Absolute ASTM D6319-00	9.45 / 240 9.1 / 230 minimum	3.2 / 0.08 1.96 / 0.05 minimum	3.94 / 0.10 1.96 / 0.05 minimum
Extra Large Aurelia Absolute ASTM D6319-00	9.45 / 240 9.1 / 230 minimum	3.2 / 0.08 1.96 / 0.05 minimum	3.94 / 0.10 1.96 / 0.05 minimum

Ordering Information

To order from your preferred distributor, specify size and item number. Packing: 200 pieces per dispenser box, 10 dispenser boxes per shipping carton.

Extra Small	Small	Medium	Large	Extra Large
98995	98996	98997	98998	98999

Authorized Aurelia Distributor:



AURELIA Div. Supermax Inc.

1899 Sequoia Drive, Aurora, IL 60506 | Toll Free: 1-877-AURELIA | Tel: 1-630-898-8886 | Fax: 1-630-898-8855
 Web: www.aureliagloves.com | Email: info@aureliagloves.com | Facebook: www.facebook.com/aureliagloves