

Perform


Aurelia® Perform™ Finger Textured
Powder Free Nitrile Examination Gloves



**Maximum Performance
Formulated for Stretch & comfort**



- Finger textured provide improved grip
- Made with soft stretch nitrile for comfort & decreased hand fatigue
- 100% visually inspected for highest quality assurance



AURELIA
Sign of Comfort™

www.aureliagloves.com | 1-877-AURELIA (1-877-287-3542)

Perform

Product Description

Aurelia® Perform™ Powder Free Nitrile Powder Free Examination Gloves are made from 100% Nitrile (Non Latex), lightweight, ambidextrous, Teal color, Finger Textured, and non-sterile. Aurelia gloves conform and are tested to the highest international standards. For more information visit www.aureliagloves.com



Aurelia® Perform™ Specifications

TENSILE TEST RESULTS (unaged)

	Tensile Strength (MPa)	Elongation at break (%)
Extra Small	25.0	700
Small	25.1	700
Medium	25.2	700
Large	25.2	700
Extra Large	25.3	7000
ASTM D6319-00	14 MPa minimum	500% minimum

AGED TENSILE TEST RESULTS (7 Days at 70° C)

	Tensile Strength (MPa)	Elongation at break (%)
Extra Small	23.0	620
Small	23.0	620
Medium	23.5	620
Large	23.5	620
Extra Large	23.5	620
ASTM D6319-00	14 MPa minimum	400% minimum

DIMENSIONS

	Length (inches/mm)	Thickness (Single-wall)	
		Palm (mils/mm)	Finger (mils/mm)
Extra Small Aurelia Perform ASTM D6319-00	9.4 / 240 8.7 / 220 minimum	2.6 / 0.066 2 / 0.05 minimum	3.94 / 0.10 2 / 0.05 minimum
Small Aurelia Perform ASTM D6319-00	9.4 / 240 8.7 / 220 minimum	2.6 / 0.066 2 / 0.05 minimum	3.94 / 0.10 2 / 0.05 minimum
Medium Aurelia Perform ASTM D6319-00	9.4 / 240 9.1 / 230 minimum	2.6 / 0.066 32 / 0.05 minimum	3.94 / 0.10 2 / 0.05 minimum
Large Aurelia Perform ASTM D6319-00	9.4 / 240 9.1 / 230 minimum	2.6 / 0.066 2 / 0.05 minimum	3.94 / 0.10 2 / 0.05 minimum
Extra Large Aurelia Perform ASTM D6319-00	9.4 / 240 9.1 / 230 minimum	2.6 / 0.066 22 / 0.05 minimum	3.94 / 0.10 2 / 0.05 minimum

ITEM NUMBER

Extra Small	Small	Medium	Large	Extra Large
92395	92396	92397	92398	92399

Ordering Information

To order from your preferred distributor, specify size and item number. Packing: 200 pieces per dispenser box, 10 dispenser boxes per shipping carton.

Authorized Aurelia Distributor:

