# Refresh

Aurelia® Refresh™ Honeycomb® Textured Peppermint Powder Free Latex Examination Gloves



Mint Scented Exclusive Honeycomb Texture



- Refreshing peppermint scent for you and your patient
- Patented Honeycomb texture for extraordinary gripping and tactile feel
- Polymer-coated enhances donning and keeps hands smooth











# Refresh

# **Product Description**

Aurelia® Refresh™ Honeycomb® Textured Powder Free Latex Examination Gloves are made from 100% Natural Malaysian rubber, ambidextrous, green color with peppermint scent, Honeycomb® Textured and non-sterile. Aurelia gloves conform and are tested to the highest international standards. For more information visit <u>www.aureliagloves.com</u>



## **Aurelia® Refresh™ Specifications**

### **TENSILE TEST RESULTS** (unaged)

### **AGED TENSILE TEST RESULTS** (7 Days at 70° C)

	Tensile Strength (MPa)	Elongation at break (%)		Tensile Strength (MPa)	Elongation at break (%)
Extra Small	30.0	970	Extra Small	28.0	900
Small	31.1	970	Small	28.0	920
Medium	31.0	970	Medium	27.5	910
Large	31.2	960	Large	27.3	890
Extra Large	31.3	970	Extra Large	27.9	900
ASTM 3578:05	18 MPa minimum	650% minimum	ASTM 3578:05	14 MPa minimum	500% minimum

### **DIMENSIONS**

	Length	Thickness (Single-wall)		
	(inches/mm)	Palm (mils/mm)	Finger (mils/mm)	
Extra Small Aurelia Refresh ASTM 3578:05	9.4 / 240 8.7 / 220 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum	
Small Aurelia Refresh ASTM 3578:05	9.4 / 240 8.7 / 220 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum	
Medium Aurelia Refresh ASTM 3578:05	9.4 / 240 9.1 / 230 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum	
Large Aurelia Refresh ASTM 3578:05	9.4 / 240 9.1 / 230 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum	
Extra Large Aurelia Refresh ASTM 3578:05	9.4 / 240 9.1 / 230 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum	

### **Ordering Information**

To order from your preferred distributor, specify size and item number. Packing: 100 pieces per dispenser box, 10 dispenser boxes per shipping carton.



Extra Small	Small	Medium	Large	Extra Large	
99225	99226	99227	99228	99229	

ITEM NUMBER

Authorized Aurelia Distributor:

CAUTION STATEMENT: This product contains natural rubber latex which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.