

Sensor Slippers™ Comfort Barriers for Digital X-Rays

Large | Yellow **REF** BWSSL / Small | Blue **REF** BWSSS


INSTRUCTIONS FOR USE

- **Large | Yellow: Fits size 2 sensors**
- **Small | Blue: Fits size 0-1 sensors**


Photos used in demonstration are Sensor Slippers™ Comfort Barriers Large, Yellow.



1
Insert digital sensor into Sensor Slipper™ Comfort Barrier.



2
Insert cushioned digital sensor into gripping-style holder of your choice.



3
Take images and discard comfort barrier after use.

CROSSTEX
A CANTHEL MEDICAL COMPANY



Crosstex International, Inc., a Cantel Medical Company
6789 W. Henrietta Road, Rush, NY 14543 USA
585.359.0130 | 631.582.6777 | crosstex.com

EC REP

MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany

CE

MD

Not made with natural rubber latex.

Made in China

BWSS-IFU-EN / REV B - 04/2021

All product names are trademarks of Crosstex International, Inc., a Cantel Medical Company, its affiliates or related companies, unless otherwise noted.

©2021 Crosstex International, Inc. All rights reserved.