



NEW



DentaMile Desk RS System

Fast, precise, and compact
chairside printing

A smile ahead together



Meet the **DentaMile Desk RS System**

Fast. Precise. Reliable.

Your fast, chairside 3D printing solution

Streamline your in-office workflow with the DentaMile Desk RS-1, designed for busy dental practices. It simplifies production, saves chair time, and reduces the need for follow-up appointments. Whether you're printing **splints, bleaching trays, crowns or surgical guides**, the DentaMile Desk RS-1 delivers exceptional accuracy at a remarkable speed, with consistent, reliable results.



TECHNICAL DATA

| | |
|-----------------------|---------------------------|
| Build space | 5.2 x 3 x 4 in. |
| Wavelength | 385 nm |
| Max. workpiece height | 90 mm |
| Resolution | Full HD, 69 µm pixel size |
| Technology | DLP with Force Feedback |
| Z resolution | 25-150 µm |
| Connectivity | WLAN, LAN, USB |
| Dimensions (W/D/H): | 12.6 x 14.8 x 15.4 in. |
| Weight | Approx. 20 lbs. |

The DentaMile Desk RS-1 is the perfect companion to the DentaMile Wash RS and DentaMile Cure RS post-processing units.



Learn more about the **DentaMile Desk RS-1** printer

FEATURES AND BENEFITS

- Fast and precise printing
- Perfect for nightguards, models, bleaching trays, crowns, and more
- Optimized for DMG LuxaPrint resins for validated results
- DMG-backed support, installation, and training
- Certified auto-calibration for consistent print quality
- Compact design ideal for chairside use
- Compatible with DentaMile connect for easy designing and direct printing
- RFID technology for automated material identification
- Open resin library