

# CaviWipes™ 1 2 3

Easy as



## CaviWipes™ HP

### Features (S.A.F.E.)

#### Simple

1-minute universal contact time. New lid design- easy to open and close securely; Dispenses wipes with precision.

#### Alcohol-free

Suitable for sensitive patient populations. Better cleaning efficacy than alcohol-based products<sup>1</sup>.

#### Fast

1-step cleaning/disinfection<sup>2</sup>. 1-handed operation lid design.

#### Effective

Protection against 67 pathogens. Fully qualified for the EPA Emerging Viral Pathogen (EVP) Claim for all virus types (enveloped, large and small non-enveloped viruses) to meet current and potentially future infection prevention needs<sup>3</sup>.

Part No.	Description
16-1100	CaviWipes HP (6" x 6.75") 160 towelettes per canister
16-1150	CaviWipes HP XL (9" x 12") 65 towelettes per canister

## CaviWipes™ 2.0

### Features

**2-Minute Universal Contact Time** for bacteria, viruses, and pathogenic fungi provided by a Quat/Low-Alcohol formulation.

**Excellent Material Compatibility** with hard, non-porous surfaces.

**1-Step Cleaning/Disinfection<sup>2</sup>** Effective against over 60 pathogens including Norovirus, SARS-CoV-2 and *Candida auris*.

#### Prepared for Today's and Tomorrow's Pathogens

Fully qualified for the EPA Emerging Viral Pathogen (EVP) Claim for all virus types (enveloped, large and small non-enveloped viruses) to meet current and potentially future infection prevention needs<sup>3</sup>.

Part No.	Description
14-1100	CaviWipes 2.0 (6" x 6.75") 160 towelettes per canister
14-1150	CaviWipes 2.0 XL (9" x 12") 65 towelettes per canister

## CaviWipes™ Bleach

### Features

**3-Minute Universal Contact Time** against 43 pathogens including *C. Diff*, TB, *C. auris*, norovirus, 8 MDROs (Multidrug Resistant Organisms).

#### No Signal Word

No Personal Protective Equipment (PPE) required, which means there is no need for gloves, mask, goggles, protective eyewear or even a face mask.

#### Minimal Residue

With our patented chemistry, our multi-purpose cleaning wipes are made of high-substrate quality and have a mild odor.

**Fully qualified for the EPA Emerging Viral Pathogen (EVP Level III) Claim** for all virus types (enveloped, large and small non-enveloped viruses) to meet current and potentially future infection prevention needs<sup>3</sup>.

Part No.	Description
13-9100	CaviWipes Bleach (6"x10.5") 90 towelettes per canister
13-1175	CaviWipes Bracket Fits all 3 CaviWipes canisters

1. Data on File.

2. Always follow product label for special pre-cleaning and PPE instructions.

3. <https://www.epa.gov/pesticide-registration/emerging-viral-pathogen-guidance-and-status-antimicrobial-pesticides>. Accessed 02.23.23.