**Composite Warmer** 



















Porcelain Veneer Tray











Restoration Tray

43 Miry Brook Road • Danbury, CT 06810 • (203) 778-0200



## **CALSET**

The Calset unit is designed to heat composite material, anesthetic solution and porcelain veneers. The Calset unit accepts composite from compules or syringed material. A tray placed on the top of the heater will either hold compules, or will have a holder for your composite dispenser. Available in all worldwide voltage.

#### **FEATURES**

- Heats composite material to 98°F (37°C), 130°F (54°C), or 155°F (68°C)
- Maintains constant temperature
- · Increases degree of cure
- Warms finishing instruments

#### **BENEFITS**

- Shortens curing time by over 80%
- Improves composite flow by 68%
- Reduces microleakage
- Decreases shrinkage stress
- Warms finishing instruments for improved marginal adaptation

There are more than 60 research studies done by leading universities worldwide on our website.

# CALSET COMPOSITE WARMER COMBINATIONS

### **CALSET WITH MULTI PRO TRAY**

- Warms a mix of composite dispensers and or syringe composite. Holds 4 compules, 2 finishing instruments.
- Includes heater base, tray and power supply

## **CALSET WITH 7 HOLE SYRINGE COMPOSITE TRAY**

- Warms 4 small syringe, and 3 large syringe
- Includes heater base, 7 Hole syringe tray and power supply.

### **CALSET WITH DISPENSER TRAY**

- Accepts composite dispenser and compules
- Includes heater base, power supply and dispenser tray with lid.

#### **CALSET WITH STANDARD TRAY**

- Warms 8 compules
- Includes heater base, power supply, standard tray and lid.

# **CALSET WITH RESTORATION TRAY**

- Holds 6 Onlays or Inlays and 2 compules
- Includes heater base, power supply, restoration tray and lid.

# **CALSET WITH VENEER TRAY**

- Warms 10 veneers
- Includes heater base, power supply, veneer tray and lid.

#### **CALSET WITH ANESTHETIC TRAY**

- Warms 3 carpules
- Includes heater base, power supply and anesthetic tray.

# TRAY ACCESSORIES AVAILABLE

- Anesthetic Tray
- 7 Hole Syringe Tray
- Veneer Tray
- Dispenser Tray
- Standard Tray
- Restoration Tray
- Multi PRO Tray





Made in the USA I www.addent.com



