

Prepare for

IMPACT

 **IMPACT AIR** *45*®



Specifically Engineered for Surgical Procedures

- + Apicoectomy (root canal drainage)
- + Odontoplasty
- + Osteoplasty (bony contouring)
- + Root Resection

Specifically Engineered for Surgical Procedures

IMPACT AIR 45® has a unique 45° angled head, allowing for maximum access and visibility. It is designed for use in surgical procedures including 3rd molar extractions, Apicoectomy (root canal drainage), Osteoplasty (bony contouring), Odontoplasty and Root Resection. It is available with fiberoptics, is fully autoclavable, has an office-replaceable turbine (minimizing downtime) and uses standard friction grip burs (featuring a push-button bur release.)

Hosing Connections (with adapter)



I.S.O. Type C
(6 Pin)



I.S.O. Type B
(4 Hole)




I.S.O. Type B
(5 Hole)

A water or saline line inside the handpiece handle directs a pure water jet (not a mist) directly onto the cutting area. Air is exhausted through the back of the handpiece.

Check out our **IMPACT AIR 45®** Accessory Line

From surgical burs to rear coupling gaskets, we offer a full line of accessories to enhance your IMPACT AIR 45® experience.



"I have been an Oral and Maxillofacial surgeon since 1985, and was trained to use straight surgical Hall drills for sectioning and removing teeth. With the development of the IMPACT AIR 45® high speed handpieces, it was a game changer in my ability to do complex surgeries with ease and speed. Having irrigation and later fiber optic lighting, I question why anyone would use anything else. It has amazing torque and reliability, and ease of repair. I unconditionally recommend it for all surgical needs."

- Dr. Ronald P Morris, DDS

PRODUCT	SKU
#401 IMPACT AIR 45® Highspeed Fiberoptic Handpiece (5 Hole)	DE-401-FO
<i>Includes: (1) IMPACT AIR 45® Fiberoptic Handpiece, (1) Head Cap Wrench</i>	
#403 IMPACT AIR 45® Highspeed Handpiece (4 Hole)	DE-403-4H
<i>Includes: (1) IMPACT AIR 45® Highspeed 4 Hole Handpiece, (1) Head Cap Wrench</i>	
#407 IMPACT AIR 45® Quick Disconnect Highspeed Fiberoptic Handpiece	DE-407-FO-QD
<i>Includes: (1) IMPACT AIR 45® Quick Disconnect Highspeed Fiberoptic Handpiece, (1) Head Cap Wrench</i>	
#408 IMPACT AIR 45® Quick Disconnect Highspeed Non-Optic Handpiece	DE-408-QD
<i>Includes: (1) IMPACT AIR 45® Quick Disconnect Highspeed Non-Optic Handpiece, (1) Head Cap Wrench</i>	

