

Aurelia® Luminance™, this specially formulated latex Glove provides excellent gripping performance. Enhanced by triple polymer washed coating & fully textured, this glove delivers greatest comfort with non slip action. Gloves are 100% visually inspected for highest quality assurance. Packed in detail with layered packing technology for easy dispensing & wrinkle free gloves.



Product Characteristics

Category	Description
Product Name	Aurelia® Luminance™
Product Type	Latex Powder Free Fully Textured Examination Gloves
Product Code	2988
Size	29885 - Extra Small 29886 - Small 29887 - Medium 29888 - Large 29889 - Extra Large
Color	Cream White
Finishes	Fully Textured
Packing	100pcs/bx, 10bxs/cs 84 Cases / Pallet

Product Specifications

Category	Description		
Material	Natural Latex		
Thickness	Location of Measurement		Single Wall (mm)
	Finger Tip (at 15+/-3 mm from the extreme tip)		min 0.08
	Palm (at centre of palm)		min 0.05
	Cuff (at 13+/-3 mm from cuff end)		min 0.05
Dimensions	Size	Palm Width (mm)	Normal (mm)
	Extra Small	80+/-5	Min 240
	Small	85+/-5	Min 240
	Medium	95+/-5	Min 240
	Large	105+/-5	Min 240
	Extra-Large	> 110	Min 240
Tensile Test Results (unaged)	Size	Tensile Strength (MPa)	Elongation at break (%)
	Extra Small	30.0	970
	Small	31.1	970
	Medium	31.1	970
	Large	31.2	970
	Extra-Large	31.3	970
	ASTM D3578	18MPa minimum	650% minimum

