

Spiral Into Control

4+1 FREE!
One Year Window*

Pivot®*Plus* for the edge you need: Improved interproximal cleaning and splatter control plus enhanced patient comfort.

New Pivot*Plus* Disposable Prophy Angles offer improved prophy efficiency while allowing you to polish adjacent teeth and interproximal areas with greater ease and patient comfort.

- Pivot*Plus* spiral cups capture more paste, allowing easier access to those hard-to-reach areas.
- Pivot*Plus* has deeper external grooves to improve splatter control compared to other spiral cups and traditional prophy angles.
- Only Pivot*Plus* spiral cups have a smooth edge that splays easily into the sulcus, improving your patient's comfort (see examples at right).

Call your favorite dental dealer to take advantage of our risk-free Pivot*Plus* trial offer, or visit www.preventech.com to request a sample. And because all Pivot DPAs are sold at the same reasonable price, Pivot*Plus* allows you to keep your hygiene supply costs under control.

PIVOT
DISPOSABLE PROPHY ANGLE



PREVENTECH®  4330-C Matthews-Indian Trail Road • Indian Trail, NC 28079 • 800.474.8681 • 704.849.2416 • fax: 704.849.2417 • preventech.com

Absolute Satisfaction Guaranteed. *Nupro Freedom HP only accepts Nupro Freedom DPAs (made in China). Young Dental and Nupro trademarks are the properties of their respective companies
*Offer expires 12/31/21. There's no limit to the number of boxes an office may order, but at least four boxes must be purchased over a one year period starting with the date of your first purchase. This offer may not be combined with any other offers.



PIVOT[®]

DISPOSABLE PROPHY ANGLE

- Designed to fit all slow-speed Doriot (straight nose cone) style handpieces
- "Runs Like Metal" performance
- Shorter profile cup for improved posterior access
- 4 out of 5 Hygienists prefer our Latex-Free Cups
- Ribbed and webbed cup design reduces splatter
- Plus Spiral Cups for improved polishing on adjacent teeth



	Item No.	Description
SOFT CUP	1100013	PIVOT DPA Gross (144 per box)
	1100017	PIVOT DPA Economy Pack (1,200 angles)
FIRM CUP	1100023	PIVOT DPA Gross (144 per box)
	1100027	PIVOT DPA Economy Pack (1,200 angles)
NL TAPERED BRUSH	1100033	PIVOT NL DPA Gross (144 per box)
	1100032	PIVOT NL DPA Bag (36 angles)
NL SOFT CUP	1100043	PIVOT NL Non-Latex DPA Gross (144 per box)
	1100047	PIVOT NL Non-Latex DPA Economy Pack (1,200 angles)
PLUS SPIRAL SOFT CUP	1133343	PIVOT NL Non-Latex DPA Gross (144 per box) NEW
	1133347	PIVOT NL Non-Latex DPA Economy Pack (1,200 angles) NEW
NL FIRM CUP	1100053	PIVOT NL Non-Latex DPA Gross (144 per box)
	1100057	PIVOT NL Non-Latex DPA Economy Pack (1,200 angles)
PLUS SPIRAL FIRM CUP	1133353	PIVOT NL Non-Latex DPA Gross (144 per box) NEW
	1133357	PIVOT NL Non-Latex DPA Economy Pack (1,200 angles) NEW

PIVOT^{NL}

CONTRA

NON-LATEX DISPOSABLE PROPHY ANGLE

- 15° Contra Bend offers greater reach, access, visibility and comfort
- Improved ergonomics helps reduce wrist fatigue
- Real savings, same price as our popular Pivot DPA
- Smooth, quiet "Runs Like Metal" performance



NL Contra SOFT CUP	1122243	PIVOT NL Contra Non-Latex DPA Gross (144 per box)
	1122247	PIVOT NL Contra Non-Latex DPA Economy Pack (1,200 angles)
NL Contra PLUS SPIRAL SOFT CUP	1144443	PIVOT NL Contra Non-Latex DPA Gross (144 per box) NEW
	1144447	PIVOT NL Contra Non-Latex DPA Economy Pack (1,200 angles) NEW
NL Contra FIRM CUP	1122253	PIVOT NL Contra Non-Latex DPA Gross (144 per box)
	1122257	PIVOT NL Contra Non-Latex DPA Economy Pack (1,200 angles)
NL Contra PLUS SPIRAL FIRM CUP	1144453	PIVOT NL Contra Non-Latex DPA Gross (144 per box) NEW
	1144457	PIVOT NL Contra Non-Latex DPA Economy Pack (1,200 angles) NEW

Buy 4 Boxes (144) or 1 Economy Pack (1200) of PIVOT, Get 1 Box (144) FREE

Name of Practice _____ Address _____ City _____ State _____ Zip _____ e-mail _____

PivotNL Contra

Specify Cup:

- Soft 1122243
 Firm 1122253

PivotNL

Specify Cup:

- Soft 1100043
 Firm 1100053

Pivot

Specify Cup:

- Soft 1100013
 Firm 1100023

PivotPLUS Contra

Specify Cup:

- Soft 1144443
 Firm 1144453

PivotPLUS

Specify Cup:

- Soft 1133343
 Firm 1133353

Just place your order for 4 Boxes or 1 Economy Pack of Pivot DPA through the dealer of your choice. Then attach your invoice to this coupon and mail it to Preventech, so we can process your FREE Box. Allow 4 to 6 weeks for delivery.

Offer expires 12/31/22. There's no limit to the number of Boxes an office may order, but at least 4 Boxes or one Economy Pack must be purchased over one year period starting with the date of your first purchase. This offer may not be combined with any other offers.

Please contact your dealer for pricing information.

Your satisfaction is 100% guaranteed. If any Preventech product does not live up to your standards, return it for a full refund. **Customer Service: 1-800-474-8681**

