

**The alloy type you want.
The consistent mix and
optimum handling you need.**





Ensures marginal integrity,
eliminates open contacts.

Premium Dispersed Phase Amalgam Alloy

	Single Spill (400 mg)		Double Spill (600 mg)		Triple Spill (800 mg)	
	Regular	Fast Set	Regular	Fast Set	Regular	Fast Set
50 Capsule Jar	61008	61009	61012	61013	61016	61017
100 Capsule Jar	61206	61207	61208	61209	61210	61211
500 Capsule Jug	61020	61021	61024	61025	61028	61029
	Green/ Gray	Green/ White	Navy Blue/ Gray	Navy Blue/ White	Yellow/ Gray	Yellow/ White

Available in 5 oz. Powder Regular 60000 Fast Set 61001

Composition Contains: 50% Mercury by Weight
Approximate Alloy Content: Silver 53%, Tin 23%, Copper 24%, Zinc <.5%



Resists fracture & ditching.

Type-T Spherical Amalgam Alloy

	Single Spill (400 mg)	Double Spill (600 mg)	Triple Spill (800 mg)
50 Capsule Jar	61039	61040	61041
100 Capsule Jar	61047	61048	61049
500 Capsule Jug	61035	61036	61037
	Aqua/Royal Blue	Aqua/Gray	Aqua/White

Composition Contains: 42.5% Mercury by Weight
Approximate Alloy Content: Silver 56%, Tin 29%, Copper 15%



Reduces early breakage.

Royale High Copper Dispersed Phase Amalgam Alloy

	Single Spill (400 mg)	Double Spill (600 mg)	Triple Spill (800 mg)
50 Capsule Jar	65007	65008	65009
100 Capsule Jar	65103	65104	65105
500 Capsule Jug	65013	65014	65015
	Green/Green	Navy Blue/Navy Blue	Yellow/Yellow

Composition Contains: 50% Mercury by Weight
Approximate Alloy Content: Silver 44%, Tin 30%, Copper 26%

Dental Amalgam For Posterior Restoration



Manufactured by: DMG America LLC • 242 South Dean Street
• Englewood, NJ 07631 • Phone: (800) 662-6383 • www.dmg-america.com