

Dux Dental Material Safety Data Sheet

Dux Dental
600 E. Hueneme Road
Oxnard, CA, USA

Phone: 800.833.8267
Fax: 805.488.2266
www.duxdental.com

Product Gingigel
Other Name Viscous 20 % Aluminum Chloride Hemostatic Gingival Retraction System
REF 13165, 13166, 13167, 13191, 13192, 13193

MSDS Number MS114

Hazardous Ingredients

Principal Hazardous Components

None

pel / tl

NA

CAS Number

NA

Percentage

NA

Physical and Chemical Characteristics

Appearance Clear, very viscous gel with a bitterseet taste.
Boiling Point > 212°F
Melting Point NA
Vapor Pressure No Data Found
Evaporation Rate Slower than water.
Odor None
Vapor Density NA
Specific Gravity 1.25
Water Solubility Aqueous System

Fire & Explosion Data

Explosion Not considered to be an explosion hazard.
Fire Not considered to be a fire hazard.
Fire Ext. Type A B C
Flashpoint NA
Special Fire Procedures None
Unusual Hazards None Known

Reactivity Data

Stability Stable
Incompatibilities None Known
Hzd Decomp Products Fumes of hydrochloric acid.
Hazardous Polymerization Will not occur.

Dux Dental Material Safety Data Sheet

Dux Dental
600 E. Hueneme Road
Oxnard, CA, USA

Phone: 800.833.8267
Fax: 805.488.2266
www.duxdental.com

Product Gingigel
Other Name Viscous 20 % Aluminum Chloride Hemostatic Gingival Retraction System
REF 13165, 13166, 13167, 13191, 13192, 13193

MSDS Number MS114

Health Hazard Information

Ingestion	May cause irritation. Dilute with water if ingested. Consult physician.
Inhalation	NA
Skin/Eye Contact	May cause irritation. Rinse with water for 15 min, consult physician.
Toxicity Data	May cause reversible or irreversible changes to exposed tissue, not permanent injury or death; can cause considerable discomfort. Use in accordance with manufacturers instructions and sound professional
Acute Hazards	None
Chronic Hazards	None
Carcinogenicity	NA, no product ingredient listed by NTP, IARC, or OSHA.

Precautions for Safe Handling and Use

Cleanup	Collect while minimizing production of dust.
Disposal	Dispose of in accordance with local, state, and federal regulations.
Precautions	As part of good industrial, personal hygiene, and safety practice avoid all unnecessary exposure to the substance and ensure prompt removal from skin, eyes, and clothing.

Occupational Control Measures

Eye Protection	NA
Gloves	NA
Dust/Particulate Mask	NA
Other	NA

Additional Information

DOT Hazard Class	Not Regulated
DOT Small Qty Exception	NA
DOT Other Exception	NA
DOT Exception Permitted	NA
Allergenic Ingredients	This product does not contain Latex, Peanut Products, Gluten, or Red Dye #3.

The information contained herein is provided in good faith and is believed to be correct as of the date hereof. However, the company makes no representation as to the comprehensiveness or accuracy of the information. It is expected that individuals receiving the information will exercise their independent judgement in determining its appropriateness for a particular purpose. Accordingly, we will not be responsible for damages of any kind resulting from the use of or reliance upon such information. No representations or warranties, either express or implied, of merchantability, fitness for a particular purpose or of any other nature are made hereunder with respect to the information set forth herein or to the product to which the information refers.