

MICRO-TOUCH[®] NITRAFREE[™] Blue

High Touch

Provides comfortable and reliable protection with low dermatitis potential

The new MICRO-TOUCH[®] NitraFree[™] Blue nitrile exam gloves provide comfortable and reliable protection with a low dermatitis potential. A thin and lightweight formulation made without sulfur-based chemical accelerants for tactility and to prevent instance of Type IV allergic reactions. The unique, berry blue coloring provides contrast while hiding blood and stains.

FEATURES & BENEFITS:

- Manufactured without sulfur-based chemical accelerators, helping to protect users from Type I latex allergies and Type IV chemical allergies
- Thin and light weight formulation provides enhanced tactile sensitivity and comfort for longer wear time
- Unique berry blue color provides contrast, while masking blood and other stains
- Tested for use with chemotherapy drugs

Recommended For

- Standard patient examinations
- Protection from latex Type I allergy and minimize chemical Type IV allergies in HCW's or patients
- Tested for use with chemotherapy drugs



MICRO-TOUCH[®] Everyday, Everywhere Protection

MICRO-TOUCH[®] premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

MICRO-TOUCH[®] NITRAFREE[™] Blue

High Touch

PRODUCT INFORMATION

Material	Nitrile	Not Made From Natural Rubber Latex	Yes
Color	Blue	Cuff Length	Standard
Powder Content	Powder-Free	External Glove Surface	Textured Fingertip
Freedom from Holes (Inspection level I)	1.5 AQL	Palm Thickness (mm/mil)	0.07 / 2.8
Finger Thickness (mm/mil)	0.10 / 3.9	Allergy Prevention	Latex (Type I), Chemical (Type IV)
Available Sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)	Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 and US FDA cleared
Silicone Free	Not Tested	Sterile	No
Antistatic	No	Double Gloving Recommendation	Yes, as underglove
Glove Length (mm/inches)	240 / 9.5	Product Segmentation	High Touch

ORDERING INFORMATION

Size	XS (5.5 - 6)	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)
Product Code	313028060	313028070	313028080	313028090	313028100

PACKAGING AND STORAGE

Packaging	200 gloves per box; 10 boxes/2000 gloves per case; *XL: 180 gloves per box; 10 boxes/1800 gloves per case
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

STANDARDS AND CERTIFICATIONS

ASTM D6319, EN 455 1-4, FDA21 CFR 177-2600-US Food Contact Approved, ISO 13485, ISO 9001

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 1-844-494-7854
Email: insidesalesus@ansell.com

Latin America & Caribbean
Tel: +52(442) 296 20 50
Email: cslac@ansell.com
Brazil Tel: +55-11-3356-3100
Brazil Email: luvas.medicas@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu

Asia Pacific
Tel: +603 8310 6688
Email: apac.medical@ansell.com

Australia & New Zealand
Tel: +61 3 9270 7270
Email: protection@ap.ansell.com