

# SciCan HYDRIM C61W G4

Instrument Washer



## Automated Instrument Washing the first step in the sterilization process



### Efficient

The time saved by automated washing compared to manual cleaning is approximately 1 hour of labor for every nine procedural set-ups<sup>1</sup>. HYDRIM G4 cycle times are as fast as 40 minutes from start to dry. HYDRIM® ensures that every set of instruments receives the same validated cleaning process.



### Capacity

HYDRIM C61W G4 has a 6 full-size cassette capacity, allowing offices to turnover 120 instruments in as little as 40 minutes, start to dry.



### Remote Data Storage & Equipment Support

Large touch screen provides in-depth cycle info and current operation. Set-up automatic send cycle data to your smart devices, computer or patient management system for easy and efficient data storage.

Allow technicians to remotely check equipment for faster, better service without having to wait for a visit.



<sup>1</sup> Tholen, Mark. "Increase production and decrease expenses in sterilization." Dental Economics, Volume 99, Issue 10, October 1, 2009.

## HYDRIM C61W G4 Specifications

### Unit Size

Dimensions (W x H x D):	598 x 520 x 526mm / 23.5" x 20.5" x 20.7"
Depth with Door Open:	82.9 cm / 32.6"
Weight:	44 kg / 97 lb

### Cycle Times\*

Rinse and Hold	5 min
Regular Wash + Dry	26 min + 10 min
Heavy Duty Wash + Dry	30 min + 10 min

### Water Consumption Per Cycle

Rinse and Hold	7.6 L / 1.7 Gal
Regular Wash	11.5 L / 3 Gal
Heavy Duty Wash	11.5 L / 3 Gal


### Electrical Connection

North America	208-240 V, 12 A, 60 Hz
Running Noise	65 dB

### Water and Drain Connections

Hot and Cold Water Tap Connection	3/4" Garden Hose, Male
Recommended Hot Water Temperature	60° C / 140° F
Cold Water Temperature	< 30° C / 86° F
Inlet Water Pressure	1-10 bar / 14.5-145 psi
Drain	2 cm / 3/4"

\*Actual cycle times will vary with incoming water temperature and pressure.



VistaPure Water Filtration shown above.

## VistaPure Water Filtration

Provides two grades of high-quality water for autoclaves, dental bottles, ultrasonic cleaners and the rinse cycle of instrument washers.

**V3000**

## HYDRIM VistaPure Connection

Protect expensive instruments and equipment by using the VistaPure system's hyperfiltered water for the final-rinse cycle of the HYDRIM for a spot-free finish.

**V3010**

## Recommended Accessories

Part Number	Description
01-113251	6-Cassette Rack
01-113252	4 Cassette Rack
01-113253	3-Cassette, 2-Basket Rack
01-113254	10-Tray, 2-Basket Rack
01-113257	Trolley Center Support
01-113545	Large Basket
01-113546	Hygiene Basket
01-113547	Basket with Hinged Lid
CS-HIPC-U	HIP ULTRA Cleaning Solution



**HYDRIM L110W G4**  
also available

### COLTENE International Dental Group

Manufactured for:

**Dent4You AG**  
Bahnhofstrasse 2  
CH-9435 Heerbrugg

Distributed by:

**SciCan Ltd.**  
1440 Don Mills Rd.  
Toronto, ON  
M3B 3P9 Canada  
T +1 416-445-1600  
customerservice@scican.com

Distributed by:

**Coltène/Whaledent Inc.**  
235 Ascot Parkway  
Cuyahoga Falls, OH  
44223 US  
T +1 330-916-8800  
service@coltene.com

VistaPure is manufactured by:

**Crosstex International, Inc.**  
1733 W. Parkside Lane  
Phoenix, AZ 85027 / USA  
T +1 800-531-2221  
VistaResearchGroup.com

[www.scican.com](http://www.scican.com)

HYDRIM is a registered trademark of SciCan Ltd. used under license by Dent4You AG. All Vista product names are trademarks of HuFriedy Mfg. Co. LLC, its affiliates or related companies, unless otherwise noted. Marks not registered in all jurisdictions.

**COLTENE**