PENTRON

# **SAFETY DATA SHEET**

Simile™ Universal Composite

Section 1. Identification				
GHS product identifier	: Simile™ Universal Composite			
Other means of identification	: Not available.			
Product type	: Paste.			
Relevant identified uses of t	the substance or mixture and uses advised against			
Product use	: Dental product: Composite			
Area of application	: Professional applications.			
Manufacturer	: Pentron Clinical 1717 West Collins Avenue Orange, CA 92867-5422 Telephone no.: 1-203-265-7397, Toll Free: 1-800-551-0283			
e-mail address of person responsible for this SDS	: edwin.varela@kavokerrgroup.com			
Emergency telephone number (with hours of operation)	: CHEMTREC® (24 hours) U.S. : 1-800-424-9300 International: +1-703-527-3887			

# Section 2. Hazards identification

OSHA/HCS status	: This material is considered hazardous by the OSHA Hazard Communication Standard (29 CFR 1910.1200).			
	Health effects are based on the uncured material.			
Classification of the substance or mixture	: SKIN IRRITATION - Category 2 EYE IRRITATION - Category 2A			
	Percentage of the mixture consisting of ingredient(s) of unknown toxicity: 96%			
GHS label elements				
Hazard pictograms				
Signal word	: Warning			
Hazard statements	: Causes serious eye irritation. Causes skin irritation.			
Precautionary statements				
Prevention	: Wear protective gloves. Wear eye or face protection. Wash hands thoroughly after handling.			
Response	: IF ON SKIN: Wash with plenty of soap and water. Take off contaminated clothing. Wash contaminated clothing before reuse. If skin irritation occurs: Get medical attention. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical attention.			
Date of issue/Date of revision	: 01/08/2015 Date of previous issue : No previous validation Version : 1 1/10			

## Section 2. Hazards identification

Storage	: Not applicable.
Disposal	: Not applicable.
Hazards not otherwise classified	: None known.

## Section 3. Composition/information on ingredients

Substance/mixture	: Mixture
Other means of identification	: Not available.

### CAS number/other identifiers

CAS number: Not applicableProduct code: Not available.			
Ingredient name	Other names	%	CAS number
1,6-hexanediyl bismethacrylate	1,6-hexanediyl bismethacrylate	5-10	6606-59-3
7,7,9(or 7,9,9)-trimethyl-4,13-dioxo-3,14-dioxa- 5,12-diazahexadecane-1,16-diyl bismethacrylate	7,7,9(or 7,9,9)-trimethyl-4, 13-dioxo-3,14-dioxa-5, 12-diazahexadecane-1, 16-diyl bismethacrylate	5-10	72869-86-4

Any concentration shown as a range is to protect confidentiality or is due to batch variation.

There are no additional ingredients present which, within the current knowledge of the supplier and in the concentrations applicable, are classified as hazardous to health and hence require reporting in this section.

## Section 4. First aid measures

Description of necessary first aid measures					
Eye contact	ntact : No special measures are required. In case of contact with eyes, rinse immediately with plenty of water. Get medical attention if symptoms occur.				
Inhalation	Inhalation : No special measures required. If inhaled, remove to fresh air. Get medical attention if symptoms occur.				
Skin contact : No special measures required. In case of contact, immediately flush skin with plenty water. Get medical attention if symptoms occur.					
Ingestion	: Wash out mouth with water. If material has been swallowed and the exposed person is conscious, give small quantities of water to drink. Stop if the exposed person feels sick as vomiting may be dangerous. Get medical attention if adverse health effects persist or are severe.				
Most important sympto	ms/effects, acute and delayed				
Potential acute health	<u>effects</u>				
Eye contact	: Causes serious eye irritation.				
Inhalation	: Exposure to decomposition products may cause a health hazard. Serious effects may be delayed following exposure.				
Skin contact	: Causes skin irritation.				
Ingestion	: No known significant effects or critical hazards.				
Over-exposure signs/s	Over-exposure signs/symptoms				
Eye contact	: Adverse symptoms may include the following:				

Inhalation	pain or irritation watering redness : No specific data.	J		
Date of issue/Date of revision	: 01/08/2015 Date of previous issue	: No previous validation	Version :1	2/10

# Section 4. First aid measures

Skin contact	: Adverse symptoms may include the following: irritation redness
Ingestion	: No specific data.
Indication of immedia	te medical attention and special treatment needed, if necessary

Notes to physician	: In case of inhalation of decomposition products in a fire, symptoms may be delayed. The exposed person may need to be kept under medical surveillance for 48 hours.
Specific treatments	: No specific treatment.
Protection of first-aiders	<ul> <li>In case of major fire and large quantities: No action shall be taken involving any personal risk or without suitable training.</li> </ul>

See toxicological information (Section 11)

## Section 5. Fire-fighting measures

	_
Extinguishing media	
Suitable extinguishing media	: Use an extinguishing agent suitable for the surrounding fire.
Unsuitable extinguishing media	: Do not use water jet.
Specific hazards arising from the chemical	: No specific fire or explosion hazard.
Hazardous thermal decomposition products	: Decomposition products may include the following materials: carbon dioxide carbon monoxide nitrogen oxides metal oxide/oxides
Special protective actions for fire-fighters	<ul> <li>In case of major fire and large quantities: Promptly isolate the scene by removing all persons from the vicinity of the incident if there is a fire. No action shall be taken involving any personal risk or without suitable training.</li> </ul>
Special protective equipment for fire-fighters	: Fire-fighters should wear appropriate protective equipment and self-contained breathing apparatus (SCBA) with a full face-piece operated in positive pressure mode.

## Section 6. Accidental release measures

Personal precautions, protec	tive equipment and emergency procedures
For non-emergency personnel	: Low release. For professional use only. Handling of product in very small amounts or in situations where release is highly unlikely
For emergency responders	: Low release. See also the information in "For non-emergency personnel".
Environmental precautions	: Low release. Avoid dispersal of spilled material and runoff and contact with soil, waterways, drains and sewers. Inform the relevant authorities if the product has caused environmental pollution (sewers, waterways, soil or air).
Methods and materials for co	ontainment and cleaning up
Small spill	: Small Quantity. For professional use only. Absorb with an inert material and place in an appropriate waste disposal container.
Large spill	: Small Quantity. For professional use only. Absorb with an inert material and place in an appropriate waste disposal container.

Date of issue/Date of revision	: 01/08/2015	Date of previous issue	: No previous validation	Version	:1	3/10
--------------------------------	--------------	------------------------	--------------------------	---------	----	------

## Section 7. Handling and storage

Precautions for safe handling	
Protective measures	: No special measures are required for small quantities under normal and intended conditions of product use. For professional use only. Put on appropriate personal protective equipment (see Section 8). Handle with care and dispose in a safe manner.
Advice on general occupational hygiene	: Eating, drinking and smoking should be prohibited in areas where this material is handled, stored and processed. Workers should wash hands and face before eating, drinking and smoking. Remove contaminated clothing and protective equipment before entering eating areas. See also Section 8 for additional information on hygiene measures.
Conditions for safe storage, including any incompatibilities	: Store in accordance with local regulations. Store in original container protected from direct sunlight in a dry, cool and well-ventilated area, away from incompatible materials (see Section 10) and food and drink. Keep container tightly closed and sealed until ready for use. Containers that have been opened must be carefully resealed and kept upright to prevent leakage. Do not store in unlabeled containers. Use appropriate containment to avoid environmental contamination.

## Section 8. Exposure controls/personal protection

Control parameters		
Occupational exposure li		
None.		
Appropriate engineering controls	No special measures are required for small quantities under normal and ir conditions of product use.	Itended
Environmental exposure controls	<ul> <li>No special measures are required for small quantities under normal and ir conditions of product use.</li> </ul>	itended
Individual protection meas	2	
Hygiene measures	<ul> <li>No special measures are required for small quantities under normal and ir conditions of product use.</li> </ul>	itended
Eye/face protection	Safety eyewear complying with an approved standard should be used whe assessment indicates this is necessary to avoid exposure to liquid splashe gases or dusts. If contact is possible, the following protection should be w the assessment indicates a higher degree of protection: chemical splash	es, mists, /orn, unless
Skin protection		
Hand protection	Chemical-resistant, impervious gloves complying with an approved standar worn at all times when handling chemical products if a risk assessment in necessary. Considering the parameters specified by the glove manufactu during use that the gloves are still retaining their protective properties. It is noted that the time to breakthrough for any glove material may be differen glove manufacturers. In the case of mixtures, consisting of several substa protection time of the gloves cannot be accurately estimated.	dicates this is rer, check should be t for different
Body protection	<ul> <li>No special measures are required for small quantities under normal and ir conditions of product use.</li> </ul>	Itended
Other skin protection	<ul> <li>Appropriate footwear and any additional skin protection measures should based on the task being performed and the risks involved and should be a specialist before handling this product.</li> </ul>	
Respiratory protection	<ul> <li>No special measures are required for small quantities under normal and ir conditions of product use.</li> </ul>	itended

# Section 9. Physical and chemical properties

### Appearance

Appearance		
Physical state	:	Solid. [Paste.]
Color	:	Various
Odor	:	Fruity ester-like
Odor threshold	:	Not available.
рН	:	Not available.
Melting point	:	Not available.
Boiling point	:	Not available.
Flash point	:	Not available.
Evaporation rate	:	Not available.
Flammability (solid, gas)	1	Not available.
Lower and upper explosive (flammable) limits	:	Not available.
Vapor pressure	1	Not available.
Vapor density	1	Not available.
Relative density	:	Not available.
Solubility	1	Insoluble in the following materials: cold water and hot water.
Solubility in water	1	Not available.
Partition coefficient: n- octanol/water	:	Not available.
Auto-ignition temperature	1	Not available.
Decomposition temperature	:	Not available.
SADT	:	Not available.
Viscosity	:	Not available.
Density	:	1.8 to 1.9 g/cm <sup>3</sup>

## Section 10. Stability and reactivity

Reactivity	: No specific test data related to reactivity available for this product or its ingredients.
Chemical stability	: The product is stable.
Possibility of hazardous reactions	: Under normal conditions of storage and use, hazardous reactions will not occur.
	Under normal conditions of storage and use, hazardous polymerization will not occur.
Conditions to avoid	: Keep away from heat and direct sunlight.
Incompatible materials	<ul> <li>Reactive or incompatible with the following materials: oxidizing materials and reducing materials. Peroxide. free radical initiators</li> </ul>
Hazardous decomposition products	: Under normal conditions of storage and use, hazardous decomposition products should not be produced.

	Date of issue/Date of revision	:01/08/2015	Date of previous issue	: No previous validation	Version	:1	5/10
--	--------------------------------	-------------	------------------------	--------------------------	---------	----	------

## Section 11. Toxicological information

#### Information on toxicological effects

Acute toxicity

Not available.

Conclusion/Summary	<ul> <li>Based on the criteria of the protocol, this product is considered non-cytotoxic per ISO 10993-5.</li> <li>Based on analysis and test results, this product is considered as biocompatible per EN ISO 7405:2008 and EN ISO 10993-1:2009.</li> </ul>
Irritation/Corrosion	

Not available.

#### **Sensitization**

Not available.

#### **Conclusion/Summary**

Skin

#### : Kligman score: Grade I (weak sensitizer)

#### **Mutagenicity**

Not available.

#### **Conclusion/Summary** : Not mutagenic in Ames test.

#### **Carcinogenicity**

Not available.

#### **Reproductive toxicity**

Not available.

#### **Teratogenicity**

Not available.

#### Specific target organ toxicity (single exposure)

Name		Route of exposure	Target organs
1,6-hexanediyl bismethacrylate	Category 3	Not applicable.	Respiratory tract irritation
7,7,9(or 7,9,9)-trimethyl-4,13-dioxo-3,14-dioxa-5, 12-diazahexadecane-1,16-diyl bismethacrylate	Category 3	Not applicable.	Respiratory tract irritation

#### Specific target organ toxicity (repeated exposure)

Not available.

#### **Aspiration hazard**

Not available.

# Information on the likely : Routes of entry anticipated: Oral, Dermal, Inhalation. routes of exposure

Potential acute health effects	
Eye contact	: Causes serious eye irritation.
Inhalation	: Exposure to decomposition products may cause a health hazard. Serious effects may be delayed following exposure.
Skin contact	: Causes skin irritation.
Ingestion	: No known significant effects or critical hazards.

#### Symptoms related to the physical, chemical and toxicological characteristics

Date of issue/Date of revision : 01/0	/08/2015 Date of previous issue	: No previous validation Version	:1 6/10
---------------------------------------	---------------------------------	----------------------------------	---------

# Section 11. Toxicological information

Eye contact	: Adverse symptoms may include the following: pain or irritation watering redness
Inhalation	: No specific data.
Skin contact	: Adverse symptoms may include the following: irritation redness
Ingestion	: No specific data.

Delayed and immediate effect	ts	and also chronic effects from short and long term exposure
Short term exposure		
Potential immediate effects	1	Not available.
Potential delayed effects	:	Not available.
Long term exposure		
Potential immediate effects	1	Not available.
Potential delayed effects	:	Not available.
Potential chronic health effe	ect	<u>s</u>
Not available.		
General	:	No known significant effects or critical hazards.
Carcinogenicity	:	No known significant effects or critical hazards.
Mutagenicity	:	No known significant effects or critical hazards.
Teratogenicity	:	No known significant effects or critical hazards.
<b>Developmental effects</b>	1	No known significant effects or critical hazards.
Fertility effects	1	No known significant effects or critical hazards.

#### Numerical measures of toxicity

Acute toxicity estimates
--------------------------

Route	ATE value
Oral	3928.5 mg/kg

## Section 12. Ecological information

#### **Toxicity**

Not available.

#### Persistence and degradability

Not available.

#### **Bioaccumulative potential**

7/10

## Section 12. Ecological information

Product/ingredient name	LogPow	BCF	Potential
7,7,9(or 7,9,9)-trimethyl-4, 13-dioxo-3,14-dioxa-5, 12-diazahexadecane-1, 16-diyl bismethacrylate	3	-	low

#### Mobility in soil

Soil/water partition: Not available.coefficient (Koc)

**Other adverse effects** : No known significant effects or critical hazards.

### Section 13. Disposal considerations

```
Disposal methods
```

```
.
```

: The generation of waste should be avoided or minimized wherever possible. Disposal of this product, solutions and any by-products should at all times comply with the requirements of environmental protection and waste disposal legislation and any regional local authority requirements.

### Section 14. Transport information

	DOT Classification	IMDG	ΙΑΤΑ
UN number	Not regulated.	Not regulated.	Not regulated.
UN proper shipping name	-	-	-
Transport hazard class(es)	-	-	-
Packing group	-	-	
Environmental hazards	No.	No.	No.
Additional information	-	-	-

Special precautions for user : Transport within user's premises: always transport in closed containers that are upright and secure. Ensure that persons transporting the product know what to do in the event of an accident or spillage.

Transport in bulk according : Not available. to Annex II of MARPOL 73/78 and the IBC Code

## Section 15. Regulatory information

I.S. Federal regulations	:	TSCA 8(a) PAIR:	mequinol				
		United States inv	ventory (TSC	A 8b): All con	nponents are I	isted or exemp	ted.
Clean Air Act Section 112 (b) Hazardous Air Pollutants (HAPs)	:	Listed					
Clean Air Act Section 602 Class I Substances	:	Not listed					
Clean Air Act Section 602 Class II Substances	:	Not listed					
DEA List I Chemicals (Precursor Chemicals)	:	Not listed					
DEA List II Chemicals (Essential Chemicals)	:	Not listed					
<u>SARA 302/304</u>							
Composition/information	on	<u>ingredients</u>					
No products were found.							
SARA 304 RQ		Not applicable.					
<u>SARA 311/312</u>							
Classification	:	Immediate (acute)	) health haza	rd			
Composition/information	on	<u>ingredients</u>					
Name		%	Fire hazard	Sudden release of pressure	Reactive	Immediate (acute) health hazard	Delayed (chronic) health hazard
1,6-hexanediyl bismethacry	late	e 5-10	No.	No.	Yes.	Yes.	No.

No.

No.

Yes.

Yes.

No.

16-diyl bismethacrylate SARA 313

Not applicable.

#### State regulations

Massachusetts	: The following components are listed: MINERAL WOOL FIBER
New York	: None of the components are listed.

5-10

- New York New Jersey
- : None of the c
  - : None of the components are listed.

### Pennsylvania

: None of the components are listed.

#### California Prop. 65


None of the components are listed.

7,7,9(or 7,9,9)-trimethyl-4,13-dioxo-3,

14-dioxa-5,12-diazahexadecane-1,

### Section 16. Other information

#### Hazardous Material Information System (U.S.A.)



Date of issue/Date of revision

## Section 16. Other information

Caution: HMIS® ratings are based on a 0-4 rating scale, with 0 representing minimal hazards or risks, and 4 representing significant hazards or risks Although HMIS® ratings are not required on SDSs under 29 CFR 1910. 1200, the preparer may choose to provide them. HMIS® ratings are to be used with a fully implemented HMIS® program. HMIS® is a registered mark of the National Paint & Coatings Association (NPCA). HMIS® materials may be purchased exclusively from J. J. Keller (800) 327-6868.

The customer is responsible for determining the PPE code for this material.

#### National Fire Protection Association (U.S.A.)



Reprinted with permission from NFPA 704-2001, Identification of the Hazards of Materials for Emergency Response Copyright ©1997, National Fire Protection Association, Quincy, MA 02269. This reprinted material is not the complete and official position of the National Fire Protection Association, on the referenced subject which is represented only by the standard in its entirety.

Copyright ©2001, National Fire Protection Association, Quincy, MA 02269. This warning system is intended to be interpreted and applied only by properly trained individuals to identify fire, health and reactivity hazards of chemicals. The user is referred to certain limited number of chemicals with recommended classifications in NFPA 49 and NFPA 325, which would be used as a guideline only. Whether the chemicals are classified by NFPA or not, anyone using the 704 systems to classify chemicals does so at their own risk.

<u>History</u>	
Date of issue/Date of revision	: 01/08/2015
Date of previous issue	: No previous validation
Version	: 1
Prepared by	: IHS
Key to abbreviations	: ATE = Acute Toxicity Estimate BCF = Bioconcentration Factor GHS = Globally Harmonized System of Classification and Labelling of Chemicals IATA = International Air Transport Association IBC = Internediate Bulk Container IMDG = International Maritime Dangerous Goods LogPow = logarithm of the octanol/water partition coefficient MARPOL 73/78 = International Convention for the Prevention of Pollution From Ships, 1973 as modified by the Protocol of 1978. ("Marpol" = marine pollution) UN = United Nations
References	: HCS (U.S.A.)- Hazard Communication Standard International transport regulations

Indicates information that has changed from previously issued version.

Notice to reader

To the best of our knowledge, the information contained herein is accurate. However, neither the above-named supplier, nor any of its subsidiaries, assumes any liability whatsoever for the accuracy or completeness of the information contained herein.

Final determination of suitability of any material is the sole responsibility of the user. All materials may present unknown hazards and should be used with caution. Although certain hazards are described herein, we cannot guarantee that these are the only hazards that exist.