

SAFETY DATA SHEET

Section 1. Product And Company Identification

Product Name: Temp-Bond NE Base

Product Use: Dental product: Temporary cement

Manufacturer: Kerr Corporation
1717 W. Collins Ave.
Orange, CA 92867-5422
U.S.A.

Information Phone Number: 1-800-841-1428 (Customer Service)

Chemical Emergency Phone Number (Chemical Spills, Leaks, Fire, Exposure or Accident only):
CHEMTREC 1-800-424-9300 (in the US) 1-703-527-3887 (Outside the US)

SDS Date of Preparation/Revision: February 20, 2019

Section 2. Hazards Identification

GHS Classification:

Aquatic Acute Toxicity Category 1

Aquatic Chronic Toxicity Category 1

Label Elements:

Warning!



Hazard Phrases

Very toxic to aquatic life with long lasting effects.

Precautionary Phrases:

Avoid release to the environment.

Collect spillage.

Dispose of contents and container in accordance with local and national regulations.

Section 3. Composition/Information on Ingredients

Component	CAS No.	Amount
Zinc oxide	1314-13-2	60-100%
White mineral oil (petroleum)	8042-47-5	5-10%

Section 4. First Aid Measures

Inhalation: Remove person to fresh air and keep comfortable for breathing. Get medical attention if you feel unwell.

Skin Contact: Gently wash with plenty of soap and water. If skin irritation occurs: Get medical attention.

Eye Contact: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Get medical attention if you feel unwell.

Ingestion: If swallowed, rinse mouth with water (only if the person is conscious). Call a POISON CENTER/doctor if you feel unwell.

Most important symptoms and effects, acute and delayed: In all cases of doubt, or when symptoms persist, seek medical attention.

Indication of immediate medical attention and special treatment, if needed: Immediate medical attention is not required.

Section 5. Fire Fighting Measures

Suitable (and Unsuitable) Extinguishing Media: Use any media appropriate for the surrounding fire. Cool fire exposed containers with water.

Specific Hazards Arising from the Chemical: Combustion may produce carbon oxides and metallic oxide.

Special Protective Equipment and Precautions for Fire-fighters: Firefighters should wear positive pressure self-contained breathing apparatus and full protective clothing for fires in areas where chemicals are used or stored. Cool fire-exposed containers with water. Contain water used in firefighting from entering sewers or natural waterways.

Section 6: Accidental Release Measures

Personal precautions, Protective equipment, and Emergency procedures: Evacuate spill area and keep unprotected personnel away. Avoid contact with eyes, skin and clothing. Wear appropriate protective clothing and equipment.

Environmental Precautions: Avoid releases to the environment. Discharging into rivers and drains is forbidden. Prevent entry to sewers and public waters. Report spill as required by local and federal regulations.

Methods and Materials for Containment and Cleaning up: Prompt cleanup and removal are necessary. Sweep and shovel into suitable containers. Absorb spills with an inert material and place in an appropriate waste disposal container. Store away from other materials.

Section 7. Handling and Storage

Precautions for Safe Handling: Prevent contact with eyes, skin and clothing. Always wear impervious gloves, chemical safety goggles and protective clothing when handling this material. Wash thoroughly

with soap and water after handling. Do not eat, drink or smoke in the work area. Use with adequate ventilation. Remove and wash contaminated clothing before reuse.

Empty containers retain product residues which can be hazardous. Follow all SDS precautions when handling empty containers.

Conditions for Safe Storage, Including any Incompatibilities: Store in a cool, dry, well-ventilated area away from direct sunlight. Keep in the original container or an approved alternative made from a compatible material, kept tightly closed when not in use. Containers that have been opened must be carefully resealed and kept upright to prevent leakage. Do not store in unlabeled containers.

Section 8. Exposure Controls / Personal Protection

Exposure Limits

Chemical	Exposure Limit
Zinc oxide	2 mg/m ³ TWA ACGIH TLV
White mineral oil (petroleum)	5 mg/ m ³ TWA OSHA PEL

Appropriate Engineering Controls: Use with adequate general or local exhaust ventilation to maintain exposure levels below the occupational exposure limits.

Respiratory Protection: None under normal use conditions with adequate ventilation. For operations where the occupational exposure limits are exceeded, an approved respirator with particulate cartridges is recommended. Equipment selection depends on contaminant type and concentration. Select in accordance with applicable regulations and good industrial hygiene practice. For firefighting, use self-contained breathing apparatus.

Hand protection: Impervious gloves are suggested to prevent skin contact. Contact your glove supplier for selection assistance.

Eye Protection: Chemical safety goggles are recommended if contact is possible.

Skin Protection: Wear protective clothing as needed to avoid skin contact and contamination of personal clothing.

Hygiene measures: Suitable eye and skin washing facilities should be available in the work area.

Section 9. Physical and Chemical Properties

Appearance:	White paste	Odor:	Odorless
Odor Threshold:	Not available	pH:	Not available
Melting/Freezing Point:	Not available	Boiling Point/Range:	Not available
Flash Point:	Not flammable	Evaporation Rate:	Not available
Flammability: (Solid, Gas)	Not applicable	Flammability Limits:	LEL: Not applicable UEL: Not applicable
Vapor Pressure:	Not available	Vapor Density:	Not available

Relative Density:	>1	Solubilities:	Insoluble in water
Partition Coefficient: (N-Octanol/Water)	Not available	Autoignition	Not available
Decomposition	Not available	Temperature:	
Temperature:		Viscosity:	Not available

Section 10. Stability and Reactivity

Reactivity: The product is not expected to be reactive.
Chemical Stability: Stable under normal storage and handling conditions.
Possibility of Hazardous Reactions: Hazardous polymerization will not occur.
Conditions to avoid: Avoid heat and direct sunlight.
Incompatible Materials: No data available.
Hazardous decomposition products: None if stored normally.

Section 11. Toxicological Information

Potential Health Effects:

Inhalation: No data available.
Skin Contact: May cause skin dryness and irritation.
Eye Contact: Causes eye irritation.
Ingestion: May be irritating to mouth, throat and stomach.
Chronic Hazards: Prolonged or repeated contact can defat the skin and lead to irritation, cracking and/or dermatitis.

Skin Sensitization: No adverse effects expected. Components are not sensitizers.

Respiratory Sensitization: No data available. This product is not expected to cause respiratory sensitization.

Germ Cell Mutagenicity: None of the components have shown mutagenic activity in animal studies.

Carcinogen: None of the components are listed as a carcinogen or potential carcinogen by IARC, NTP, ACGIH, or OSHA.

Developmental / Reproductive Toxicity: None of the components have been shown to cause reproductive or developmental toxicity.

Specific Target Organ Toxicity (Single Exposure): No data available.

Specific Target Organ Toxicity (Repeated Exposure): No data available.

Aspiration Toxicity:

White mineral oil (petroleum) is an aspiration hazard category 1.

Acute Toxicity Values:

Zinc oxide: LD50 Oral rat: >15000 mg/kg; LD50 Dermal rat: >2000 mg/kg;
 LC50 Inhalation rat: 5.7 mg/L/4hr

White mineral oil (petroleum): LD50 Oral rat: >5000 mg/kg; LD50 Dermal rabbit: >2000 mg/kg; LC50 Inhalation rat: >5000 mg/L/4hr

Section 12. Ecological Information

Toxicity:

Zinc oxide: 96 hr LC50 Rainbow trout 1.1 mg/L; 48 hr EC50 Daphnia magna 98 µg/L;

72 hr EC50 Algae 0.042 mg/L

White mineral oil (petroleum): 96 hr LC50 Lepomis macrochirus >10000 mg/L

This product is classified as very toxic to the aquatic environment with long-term adverse effects. Releases to the environment should be avoided.

Persistence and degradability: Biodegradation is not applicable to inorganic substances.

Bioaccumulative Potential:

Zinc oxide has a BCF of 60960, potential for bioaccumulative is high.

White mineral oil (petroleum): log P_{ow} >6, potential for bioaccumulative is high.

Mobility in Soil: No data available.

Other Adverse Effects: No data available.

Section 13. Disposal Considerations

Disposal: For unused product, dispose of in accordance with Federal and local regulations.

Container Disposal: Dispose of empty container in accordance with Federal and local regulations.

Section 14. Transport Information

	UN Number	UN Proper Shipping Name	Hazard Class(s)	Packing Group	Environmental Hazards
US DOT	UN3077	Environmentally hazardous substances, solid, n.o.s. (zinc oxide). Marine pollutant (zinc oxide)	9	III	Yes
EU ADR/RID	UN3077	Environmentally hazardous substance, solid, n.o.s. (zinc oxide)	9	III	Yes
IMDG	UN3077	Environmentally hazardous substance, solid, n.o.s. (zinc oxide). Marine pollutant (zinc oxide)	9	III	Yes
IATA/ICAO	UN3077	Environmentally hazardous substance, solid, n.o.s. (zinc oxide)	9	III	Yes

Special Precautions for User: None identified

Transport in Bulk According to Annex II MARPOL 73/78 and the IBC Code: Not applicable – product is transported only in packaged form.

Section 15. Regulatory Information

U.S. Federal Regulations:

EPA SARA 311/312 Hazard Classification: Refer to Section 2 for OSHA Hazard Classification.

EPA SARA 313: This Product Contains the Following Chemicals Subject to Annual Release Reporting Requirements Under SARA Title III, Section 313 (40 CFR 372): None

Protection Of Stratospheric Ozone: This product is not known to contain or to have been manufactured with ozone depleting substances as defined in 40 CFR Part 82, Appendix A to Subpart A.

CERCLA SECTION 103: This product is not subject to CERCLA reporting requirements; however, many states have more stringent release reporting requirements. Report spills required under federal, state and local regulations.

International Inventories

US EPA TSCA Inventory: All of the components of this product are listed on the Toxic Substances Control Act (TSCA) Chemical Substances Inventory or exempt.

Canada CEPA: All of the components of this material are listed on the DSL or exempt.

Section 16. Other Information

Effective Date: February 20, 2019

Supersedes Date: February 11, 2015

Revision Summary: All Sections – New SDS format

The information and recommendations set forth herein are taken from sources believed to be accurate as of the date of preparation, however, KERR Corporation makes no warranty with respect to the accuracy or suitability of the recommendations, and assumes no liability to any use thereof.